

Test kits for water analysis

Description of individual parameters and tests

Test kits for water analysis

Contents

Test kits for water analysis	50
<i>VISOCOLOR® alpha</i>	51
<i>VISOCOLOR® ECO</i>	52
<i>VISOCOLOR® HE</i>	53
<i>VISOCOLOR®</i> program	54
Analytical principles	56
Description of individual parameters and tests.....	56
Reagent cases.....	70
Compact photometer PF-12	74

visocolor®

Test kits for water analysis

VISOCOLOR® test kits and reagent cases

Compact and flexible

- Complete mini-laboratories with reagents and accessories for water analysis
- Chemical analysis without additional accessories and without the need for any prior experience in chemistry
- Suitable for analysis in labs or directly on-site
- 3 product lines with different accuracies, precisions and sensitivities for universal use depending on the analytical requirement
- Various measuring methods and detection principles for all parameters from acidity to zinc
- VISOCOLOR® reagent cases as portable laboratories with individual combinations of different test kits

Simple and precise

- Comfortable test procedures since test kits are based on simple chemical-analytical methods like colorimetry and titration
- Instructions in different languages and with pictograms for safe and simple test performance
- Color-coded reagent bottles for clear identification of reagents
- Fast-dissolving reagents save time and facilitate the daily work – no crushing of tablets and no stirring of reaction solutions are necessary

Reliable and safe

- Reliable and comparable results – the reaction principles of VISOCOLOR® tests are based on internationally acknowledged regulations like DIN-, EN-, ISO- and EPA standards
- Maximum safety for the user and easy disposal of used reagents by avoiding of dangerous and environmentally hazardous substances
- Low susceptibility to interferences, high selectivity with regards to the substance to be analyzed as well as compensation of turbidities and colors guarantee reliable results
- Additional increase of accuracy and reproducibility by photometric determination of VISOCOLOR® ECO tests with the photometer PF-12



VISOCOLOR® alpha tests use colorimetric, as well as titrimetric procedures. Use of multicomponent reagents results in a very convenient, rapid and safe handling, because often only one reagent is needed for each test.

The reagent bottles are packed in practical blister packs. The cardboard back is used for opening and closing the package, and also provides all information required for the test: instruc-

tions for use in 6 languages with pictograms, as well as a color comparison chart for colorimetric evaluations. The blister packs of the VISOCOLOR® alpha test kits have a punched hole for convenient storage or display in showrooms or at sales counters.

Colorimetric tests

Principle:


Colorimetry with color comparison card

- Visual evaluation
- Environment-friendly
- Cost-efficient
- Convenient handling, as easy as test strips
- Accurate results
- Handy packages
- With pictogram instructions
- Reagent bottles with clear dosing instructions




Test kit consists of plastic pack with:


- Sample tube with 5 mL ring mark
- Color coded bottles with liquid or powder reagents
- Measuring spoon for accurate dosage of solid reagents
- Color scale with at least 5 gradations




Fill the sample



Add reagent



Mix



Analyze

Titration test kits

Principle:


Titration with drop counting

- Visual evaluation
- Environment-friendly
- Cost-efficient
- Convenient handling, as easy as test strips
- Accurate results
- Indicator and titration solution in one dropping bottle
- Handy packages
- With pictogram instructions
- Reagent bottles with clear dosing instructions




Test kit consists of plastic pack with:


- Sample tube with 5 mL ring mark
- One dropping bottle with mixture of indicator and titration solution



Fill the sample



Add reagent



Color change

Count the drops: 1 drop = 1 measuring unit

Test kits for water analysis

VISOCOLOR® ECO

VISOCOLOR® ECO presents a product group of colorimetric and titrimetric test kits which avoids hazardous substances wherever possible. With VISOCOLOR® ECO even water constituents with low limiting values can be determined with suffi-

cient accuracy. All VISOCOLOR® ECO test kits are packed in an environment-friendly box and contain easy to understand instructions in 6 languages.

Colorimetric test kits

Principle:


Colorimetry with color comparison card

- Visual and photometric evaluation (PF-12)
- Environment-friendly
- Economic
- Convenient handling
- Higher accuracy and sensitivity
- With pictogram instructions
- Reagent bottles with clear dosing instructions
- Compensation of turbidity and colors
- Refill packs available



Test kit consists of cardboard box with:


- 2 measuring tubes 20 mm diameter with screw caps
- Holder for the measuring tubes
- Color coded bottles with liquid or powder reagents
- Graduated plastic syringe 5 mL for convenient sample dosage
- Measuring spoon for accurate dosage of solid reagents
- Color comparison card with at least 5 gradations




Fill the sample




Add reagent



Mix



Wait



Analyze

Titration test kits

Principle:


Titration with drop counting

- Visual evaluation
- Environment-friendly
- Economic
- Convenient handling
- Higher accuracy and sensitivity
- Clearer color change due to separated dropping reagents
- Reagent bottles with clear dosing instructions




Test kit consists of cardboard box with:


- Sample tube with 5 mL ring mark
- Graduated plastic syringe 5 mL for convenient sample dosage
- Dropping bottle(s) with indicator solution
- Dropping bottle(s) with titration solution




Fill the sample



Add indicator



Titration solution



Color change

VISOCOLOR® HE test kits are highly sensitive colorimetric test kits. In comparison with conventional VISOCOLOR® kits, their sensitivity is enhanced by increasing the length of the test tube and the use of highly sensitive reagents. This technique allows a 10-fold to 100-fold increase in sensitivity. Each VISOCOLOR® HE test kit is packed in a robust plastic box

Colorimetric test kits

Principle:


High sensitivity colorimetry with comparator block and color comparison disc

- Visual evaluation
- Environment-friendly
- Convenient handling
- Highest accuracy due to extremely narrow gradation
- Highest sensitivity down to 0.002 mg/L due to longer measuring tubes
- Reagent bottles with clear dosing instructions
- Compensation of turbidity and colors
- Refill packs available




Test kit consists of plastic box with:


- 2 measuring tubes 20 mm diameter with screw caps
- Comparator block with color comparison disk
- Color coded bottles with liquid or powder reagents
- Measuring spoon for accurate dosage of solid reagents
- Beaker for convenient sample dosage




Fill the sample




Add reagent



Mix



Wait



Analyze

which contains the comparator block with color comparison disk and all required reagents. The VISOCOLOR® titration test kits are based on principles of volumetric analysis. On the graduated syringe the results can be read off in mg/L or in alternative dimensions.

Titration test kits

Principle:


High sensitive volumetric analysis with graduated syringe

- Visual evaluation
- Convenient handling
- Highest accuracy due to narrowly graduated syringe
- Reagent bottles with clear dosing instructions
- Clearer color change due to separated dropping reagents
- Refill packs available




Test kit consists of cardboard box with:


- Sample tube with 5 mL ring mark
- Graduated syringe for precise reagent dosage
- Bottle(s) with indicator solution
- Bottle(s) with titration solution




Fill the sample




Add reagent



Mix



Wait



Analyze

Test kits for water analysis

VISOCOLOR® program

VISOCOLOR® tests kits are complete plastic boxes with all reagents and accessories for a test.

Refill packs are for the replacement of used reagents in a

test kit or reagent case. They don't come with a color chart or accessories.

Ordering information

Test	Range	Type	No of tests	REF test kit	refill pack
Acidity AC 7 (base capacity)	0.2–7.0 mmol/L H ⁺ ¹⁾	HE	200	915 006	915 206
Alkalinity AL 7* (total)	0.2–7.0 mmol/L OH ⁻ ¹⁾	HE	200	915 007	915 207
Alkalinity (p/m-Wert) see Carbonate hardness C 20					
Aluminum	0.10–0.50 mg/L Al ³⁺	ECO	50	931 006	931 206
Ammonium 15*	0.5–15 mg/L NH ₄ ⁺	ECO	50	931 010	931 210
Ammonium	0.2–3 mg/L NH ₄ ⁺	alpha	50	935 012	–
Ammonium 3*	0.2–3 mg/L NH ₄ ⁺	ECO	50	931 008	931 208
Ammonium*	0.02–0.50 mg/L NH ₄ ⁺	HE	110	920 006	920 106
Calcium CA 20*	0.6–25.0 °e / 0.1–3.6 mmol/L Ca ²⁺ ¹⁾	HE	200	915 010	915 210
Calcium*	1 drop \triangleq 5 mg/L Ca ²⁺	ECO	100	931 012	–
Carbonate hardness	1 drop \triangleq 1.25 °e \triangleq 17.8 mg/L CaCO ₃	alpha	100	935 016	–
Carbonate hardness	1 drop \triangleq 1.25 °e \triangleq 17.8 mg/L CaCO ₃	ECO	100	931 014	–
Carbonate hardness C 20 (p-/m-Wert)	0.6–25.0 °e / 0.2–7.2 mmol/L H ⁺ ¹⁾	HE	200	915 003	915 203
Chloride*	1–60 mg/L Cl ⁻	ECO	90	931 018	931 218
Chloride CL 500*	5–500 mg/L Cl ⁻ ¹⁾	HE	300	915 004	915 204
Chlorine, free	0.25–2.0 mg/L Cl ₂	alpha	150	935 019	–
Chlorine 2, free + total	0.1–2.0 mg/L Cl ₂	ECO	150	931 015	931 215
free Chlorine 2	0.1–2.0 mg/L Cl ₂	ECO	150	931 016	931 216
Chlorine 6, free + total ²⁾	0.05–6.00 mg/L Cl ₂	ECO	200	–	931 217
free Chlorine 6 ²⁾	0.05–6.00 mg/L Cl ₂	ECO	400	–	931 219
Chlorine	0.02–0.60 mg/L Cl ₂	HE	160	920 015	920 115
Chlorine + pH see Swimming pool					
Chlorine dioxide* NEW!	<0.2–3.8 mg/L ClO ₂	ECO	150	931 021	931 221
Chromium(VI)*	0.02–0.50 mg/L Cr(VI)	ECO	140	931 020	931 220
Copper	0.1–1.5 mg/L Cu ²⁺	ECO	100	931 037	931 237
Copper	0.04–0.50 mg/L Cu ²⁺	HE	150	920 050	920 150
Cyanide*	0.01–0.20 mg/L CN ⁻	ECO	100	931 022	931 222
Cyanide*	0.002–0.04 mg/L CN ⁻	HE	55	920 028	920 128
Cyanuric acid	10–100 mg/L Cya	ECO	100	931 023	931 223
DEHA (diethylhydroxylamine)	0.01–0.30 mg/L DEHA	ECO	125	931 024	931 224
Fluoride ³⁾	0.1–2.0 mg/L F ⁻	ECO	150	–	931 227
total Hardness*	1 drop \triangleq 1.25 °e \triangleq 17.8 mg/L CaCO ₃	alpha	100	935 042	–
total Hardness	1 drop \triangleq 1.25 °e \triangleq 17.8 mg/L CaCO ₃	ECO	110	931 029	–
total Hardness H 20 F	0.6–25.0 °e / 0.1–3.6 mmol/L Ca ²⁺ ¹⁾	HE	200	915 005	915 205
total Hardness H 2*	0.05–2.00 °d / 0.01–0.36 mmol/L Ca ²⁺ ¹⁾	HE	200	915 002	915 202
residual Hardness *	0.04–0.30 °d	alpha	200	935 080	–
Hydrazine*	0.05–0.40 mg/L N ₂ H ₄	ECO	130	931 030	931 230
Iron 1* NEW!	0.04–1.0 mg/L Fe	ECO	200	931 025	931 225
Iron 2	0.04–1.0 mg/L Fe	ECO	100	931 026	931 226
Iron	0.01–0.20 mg/L Fe	HE	300	920 040	920 140
Manganese*	0.1–1.5 mg/L Mn	ECO	70	931 038	931 238
Manganese*	0.03–0.50 mg/L Mn	HE	100	920 055	920 155

¹⁾ For titration test kits the range can be increased with additional reagent syringe.

²⁾ only for the photometric determination with PF-11, PF-12 and VISOCOLOR® photino

³⁾ only for the photometric determination with PF-11 and PF-12

⁴⁾ based on the chemical procedures of the German Standard Methods (DEV)

* This product contains harmful substances which must be specially labeled as hazardous. For detailed information please see MSDS.

Test kits for water analysis

VISOCOLOR® program

Test	Range	Type	No of tests	REF test kit	refill pack
Nickel*	0.1–1.5 mg/L Ni ²⁺	ECO	150	931 040	931 240
Nitrate*	2–50 mg/L NO ₃ ⁻	alpha	100	935 065	–
Nitrate*	1–120 mg/L NO ₃ ⁻	ECO	110	931 041	931 241
Nitrite*	0.05–1.0 mg/L NO ₂ ⁻	alpha	200	935 066	–
Nitrite	0.02–0.5 mg/L NO ₂ ⁻	ECO	120	931 044	931 244
Nitrite	0.005–0.10 mg/L NO ₂ ⁻	HE	150	920 063	920 163
pH 5–9*	pH 5.0–9.0	alpha	200	935 075	–
pH 4.0–9.0	pH 4.0–9.0	ECO	450	931 066	931 266
pH 4.0–10.0	pH 4.0–10.0	HE	500	920 074	920 174
pH 6.0–8.2 ²⁾	pH 6.0–8.2	ECO	150	–	931 270
Phosphate*	2–20 mg/L PO ₄ ³⁻	alpha	70	935 079	–
Phosphate*	0.2–5 mg/L P	ECO	80	931 084	931 284
Phosphate*	0.05–1.0 mg/L P	HE	300	920 082	920 182
Phosphate* (DEV) ⁴⁾	0.01–0.25 mg/L P	HE	100	920 080	920 180
Potassium*	2–15 mg/L K ⁺	ECO	60	931 032	931 232
Residual hardness	see Hardness (residual)				
Oxygen*	1–10 mg/L O ₂	ECO	50	931 088	931 288
Oxygen SA 10*	0.2–10 mg/L O ₂ ¹⁾	HE	100	915 009	915 209
Silica / silicon	0.2–3.0 mg/L SiO ₂	ECO	80	931 033	931 233
Silica* / silicon	0.01–0.30 mg/L Si	HE	120	920 087	920 187
Sulfate*	25–200 mg/L SO ₄ ²⁻	ECO	100	931 092	931 292
Sulfide*	0.1–0.8 mg/L S ²⁻	ECO	90	931 094	931 294
Sulfite	1 drop \triangleq 1 mg/L SO ₃ ²⁻	ECO	60	931 095	–
Sulfite SU 100*	2–100 mg/L SO ₃ ²⁻ ¹⁾	HE	100	915 008	915 208
Swimming pool (Chlorine + pH)	0.1–2.0 mg/L Cl ₂ pH 6.9–8.2	ECO	150 150	931 090	931 290
Zinc	0.5–3 mg/L Zn ²⁺	ECO	120	931 098	931 298

¹⁾ For titration test kits the range can be increased with additional reagent syringe.

²⁾ only for the photometric determination with PF-11, PF-12 and VISOCOLOR® photino

³⁾ only for the photometric determination with PF-11 and PF-12

⁴⁾ based on the chemical procedures of the German Standard Methods (DEV)

* This product contains harmful substances which must be specially labeled as hazardous. For detailed information please see MSDS.



Test kits for water analysis

Analytical principles

Colorimetry

In colorimetric procedures advantage is taken of the fact that certain reagents form colored compounds with the substances to be determined. The intensity of the color is directly related to the concentration of the substance in question. For example in the case of the *VISOCOLOR*[®] Nitrite, the reagents form a blue-red dye with nitrite. The resulting color intensity is proportional to the concentration of nitrite. In the case of pH measurements the use of specific indicator mixtures permits the formation of a characteristic color for each pH value. The reaction colors obtained are compared with a range of standards in a sample vessel called comparator. As soon as the reaction color has been matched to a comparison color, the result can be read from the comparator or color charts.

Volumetric analysis (titration)

For certain substances it is difficult or even impossible to convert them to compounds which can be colorimetrically evaluated. In many of these cases titrimetric methods are used for analysis. This measuring principle is based on the following idea: In volumetric analysis, individual drops of a titration solution are added to a defined volume of sample solution. The active substance in the titration solution reacts with the substance to be determined in the sample. After complete reaction further addition of titration solution would cause an excess of the active substance. The point of complete reaction (end point or equivalence point) is visualized by the color change of an indicator added to the sample.

Description of individual parameters and tests

Acidity

H⁺

In natural unpolluted waters mainly carbonic acid but also humic acids are present. With this test kit, all acids, i. e. also those present in industrial waters, can be determined.

Reaction basis:

Titrimetric determination of acids with sodium hydroxide solution against the p-indicator (Reaction basis according to DIN 38409-H7).

Note:

To differentiate between mineral acids and carbonic acid, the sample should be titrated against the m-indicator from the test kit *VISOCOLOR*[®] HE Carbonate hardness C 20.

VISOCOLOR[®] HE Acidity AC 7

REF 915 006

Refill pack

REF 915 206

Type: titration test kit
Range: 0.2–7.0 mmol/L H⁺
1 gradation mark = 0.2 mmol/L H⁺
Sufficient for: about 200 tests with an average acid content of 4 mmol/L H⁺
Shelf life: at least 2 years
Sea water suitability: yes

Alkalinity (total)

OH⁻

All compounds which cause a pH increase above pH 7 are determined, e.g. hydroxide, carbonate, hydrogen carbonate etc.

Reaction basis:

Titrimetric alkalinity determination with hydrochloric acid against the m-indicator (Reaction basis according to DIN 38409-H7).

Note:

To differentiate between hydroxide, carbonate and hydrogen carbonate you should use the test kit *VISOCOLOR*[®] HE Carbonate hardness C 20 (see German Standard Methods DIN 38 409-H7).

VISOCOLOR[®] HE Alkalinity AL 7

REF 915 007

Refill pack

REF 915 207

Type: titration test kit
Range: 0.2–7.0 mmol/L OH⁻
1 gradation mark = 0.2 mmol/L OH⁻
Sufficient for: about 200 tests with an average acid content of 4 mmol/L OH⁻
Shelf life: at least 2 years
Sea water suitability: yes

Aluminum

Al³⁺

Aluminum is the most common metal in our anthroposphere and behind oxygen and silica the most common element of the earth's crust. Because of its big affinity to oxygen, aluminum does not exist in elemental form in nature, but in different oxidized compounds.

For drinking water, the WHO recommends a threshold value of 0.2 mg/L Al³⁺. In accordance to the EU council guideline 98/83/EEC, the threshold value for drinking water is 0.2 mg/L Al³⁺. In natural waters, the concentration of aluminum compounds is usually low, but waste water can contain aluminum in higher concentrations, e.g. at electroplating companies or paper mills. Different national regulations tolerate 2–3 mg/L Al³⁺ in effluents from various industries (metal, electroplating and printing industries).

Reaction basis:

Colorimetric determination with chromazurol S.

VISOCOLOR[®] ECO Aluminum

REF 931 006

Refill pack

REF 931 206

Type: colorimetric test kit
Range: 0 · 0.10 · 0.15 · 0.20 · 0.25 · 0.30 · 0.40 · 0.50 mg/L Al³⁺
Sufficient for: 50 tests
Shelf life: at least 2 years
Sea water suitability: yes, after dilution (1+9)

Test kits for water analysis

Description of individual parameters and tests

Ammonium



Ammonium ions occur primarily in domestic waste waters, but frequently also in industrial waste waters. In surface and ground waters, ammonium ions indicate decomposition of animal or vegetable matter. Control of the ammonium values is therefore important to the water supply, as a contamination indicator.

Reaction basis:

DEV procedure: in alkaline medium ammonium ions react with chlorine to form chloroamines which in the presence of a catalyst form a blue colored indophenol with phenols (Reaction basis according to DIN 38406-E5).

Primary amines react like ammonium ions resulting in high results.

Chlorine-consuming substances can lower the result or even inhibit the reaction depending on their concentration.

VISOCOLOR® alpha Ammonium

REF 935 012

Type: colorimetric test kit
Range: $0 \cdot 0.2 \cdot 0.5 \cdot 1 \cdot 2 \cdot 3 \text{ mg/L NH}_4^+$
Sufficient for: 50 tests
Shelf life: at least 1.5 years
Sea water suitability: yes, after dilution (1+9)

VISOCOLOR® ECO Ammonium 15

REF 931 010

Refill pack

REF 931 210

Type: colorimetric test kit
Range: $0 \cdot 0.5 \cdot 1 \cdot 2 \cdot 3 \cdot 5 \cdot 7 \cdot 10 \cdot 15 \text{ mg/L NH}_4^+$
Sufficient for: 50 tests
Shelf life: at least 1.5 years
Sea water suitability: yes, after dilution (1+9)



VISOCOLOR® ECO Ammonium 3

REF 931 008

Refill pack

REF 931 208

Type: colorimetric test kit
Range: $0 \cdot 0.2 \cdot 0.3 \cdot 0.5 \cdot 0.7 \cdot 1 \cdot 2 \cdot 3 \text{ mg/L NH}_4^+$
Sufficient for: 50 tests
Shelf life: at least 1.5 years
Sea water suitability: yes, after dilution (1+9)

VISOCOLOR® HE Ammonium

REF 920 006

Refill pack

REF 920 106

Type: highly sensitive test kit
Range: $0.0 \cdot 0.02 \cdot 0.04 \cdot 0.07 \cdot 0.10 \cdot 0.15 \cdot 0.20 \cdot 0.30 \cdot 0.40 \cdot 0.50 \text{ mg/L NH}_4^+$
Sufficient for: 110 tests
Shelf life: at least 1 year
Sea water suitability: no

Bromine



Bromine and brominating reagents such as 1,3-dibromo-5,5-dimethylhydantoin (DBH) are used – like chlorine – for disinfecting swimming pool water. For determination of bromine you may use all VISOCOLOR® Chlorine tests (see page 59). A factor for conversion is provided in the instructions.

Calcium



Calcium is widely distributed in nature in rocks and in water. Water containing calcium and magnesium causes problems in industry as well as in households, since during boiling calcium carbonate precipitate as detrimental boiler scale. Additionally, calcium ions inhibit soap foaming.

Reaction basis:

Complexometric titration after precipitation of the magnesium salts (Reaction basis according to DIN 38406-E3).

VISOCOLOR® ECO Calcium

REF 931 012

Type: titration test kit
Range: 1 drop = 5 mg/L Ca^{2+}
Sufficient for: about 100 tests with an average calcium concentration of 50 mg/L Ca^{2+}
Shelf life: at least 1.5 years
Sea water suitability: yes, after dilution (1+4)



VISOCOLOR® HE Calcium CA 20

REF 915 010

Refill pack

REF 915 210

Type: titration test kit
Range: $0.6\text{--}25.0 \text{ }^\circ\text{e}$ or $0.1\text{--}3.6 \text{ mmol/L Ca}^{2+}$
Sufficient for: about 200 tests with an average calcium hardness of $12.5 \text{ }^\circ\text{e}$ or $1.8 \text{ mmol/L Ca}^{2+}$
Shelf life: at least 2 years
Sea water suitability: yes, after dilution (1+4)

Test kits for water analysis

Description of individual parameters and tests

Carbonate hardness



Carbonate hardness is the part of calcium and magnesium ions which are present as carbonate or hydrogen carbonate.

Reaction basis:

The determination is preformed as titration with hydrochloric acid against a mixed indicator, which changes color at pH 4.5 (Reaction basis according to DIN EN ISO 9963-1 C24).

Normally the carbonate hardness is smaller than the total hardness. If the carbonate hardness is larger than the total hardness, you should trace the origin of this abnormal conditions, e.g. discharge of alkali hydrogen carbonates or high buffer capacity.

VISOCOLOR® alpha Carbonate hardness REF 935 016

Type: titration test kit
Range: 1 drop \triangleq 1.25 °e \triangleq 17.8 mg/L CaCO_3
Sufficient for: about 100 tests with an average hardness of 12.5 °e
Shelf life: at least 1.5 years
Sea water suitability: yes

VISOCOLOR® ECO Carbonate hardness REF 931 014

Type: titration test kit
Range: 1 drop \triangleq 1.25 °e \triangleq 17.8 mg/L CaCO_3
Sufficient for: about 100 tests with an average hardness of 12.5 °e
Shelf life: at least 2 years
Sea water suitability: yes



VISOCOLOR® HE Carbonate hardness C 20

REF 915 003
REF 915 203

Refill pack

Type: titration test kit
Range: 0.6–25.0 °e or 0.2–7.0 mmol/L H^+
Sufficient for: about 200 tests with an average hardness of 12.5 °e or 3.6 mmol/L H^+
Shelf life: at least 2 years
Sea water suitability: yes

With this test kit you can also determine the partial alkalinity (p-value) in addition to the carbonate hardness (m-value).

Carbonic acid



Carbonic acid is a natural component of water acidity. Determination is performed with the test kit VISOCOLOR® Acidity AC 7 (see page 56).

Chloride



Chloride ions occur in all natural waters. Their concentration depends on geological and regional factors. In waste waters and polluted rivers, chloride concentrations can reach high values.

Reaction bases:

- (a) Mercurimetric titration
- (b) Mercury(II) thiocyanate method

VISOCOLOR® ECO Chloride

REF 931 018

Refill pack

REF 931 218

Type: colorimetric test kit
Range: 1 · 2 · 4 · 7 · 12 · 20 · 40 · 60 mg/L Cl^-
Reaction basis: (b) Mercury(II) thiocyanate method
Sufficient for: 90 tests
Shelf life: at least 1 year
Sea water suitability: no

VISOCOLOR® HE Chloride CL 500

REF 915 004

Refill pack

REF 915 204

Type: titration test kit
Range: 5–500 mg/L Cl^-
1 gradation mark = 5 mg/L Cl^-
Reaction basis: (a) Mercurimetric titration
Sufficient for: about 300 tests with an average chloride ion concentration of 200 mg/L Cl^-
Shelf life: at least 2 years
Sea water suitability: yes, after dilution (1:50)



Test kits for water analysis

Description of individual parameters and tests

Chlorine

Cl₂

The addition of chlorine to swimming pools, water reservoirs and water mains is an approved procedure to rid the water of germs. With the correct dose harmful microorganisms are destroyed, many impurities removed and the growth of algae is prevented. However it is imperative that the chlorine content is regularly checked, since excessive chlorine not only impairs the smell and taste of the water but can be hazardous. One distinguishes between free chlorine and bound chlorine (chloroamines); the sum of both is called total chlorine.

Reaction basis:

At a pH value of 5–6 free chlorine reacts with *N,N*-diethyl-1,4-phenylene diamine (DPD) to form a red-violet dye. In the presence of iodide ions the content of total chlorine can be determined, too. (Reaction basis according to DIN ISO 7393 G4-2).

The determination of free chlorine includes the concentrations of bromine, bromoamine, chloroamine, iodine and in part chlorine dioxide.

1.0 mg/L Cl₂ \triangleq 2.3 mg/L Br₂ \triangleq 3.6 mg/L I₂

Higher-valency manganese compounds simulate free chlorine.

Note:

When measuring the chlorine content in swimming pools it is recommended to check the pH value as well. For this purpose we supply combination test kits *VISOCOLOR*[®] *ECO* Swimming pool (page 69).

VISOCOLOR[®] *alpha* Chlorine, free REF 935 019

Type: colorimetric test kit
Range: 0.25 · 0.5 · 1.0 · 1.5 · 2.0 mg/L Cl₂
Sufficient for: 150 tests
Shelf life: at least 1.5 years
Sea water suitability: yes

VISOCOLOR[®] *ECO* Chlorine 2, free + total REF 931 015 Refill pack REF 931 215

Type: colorimetric test kit
Range: < 0.1 · 0.1 · 0.2 · 0.3 · 0.4 · 0.6 · 0.9 · 1.2 · 2.0 mg/L Cl₂
Sufficient for: 150 tests
Shelf life: at least 1.5 years
Sea water suitability: yes



VISOCOLOR[®] *ECO* free Chlorine 2 Refill pack

REF 931 016
REF 931 216

Type: colorimetric test kit
Range: < 0.1 · 0.1 · 0.2 · 0.3 · 0.4 · 0.6 · 0.9 · 1.2 · 2.0 mg/L Cl₂
Sufficient for: 150 tests
Shelf life: at least 1.5 years
Sea water suitability: yes

VISOCOLOR[®] *ECO* Chlorine 6, free + total REF 931 217

Type: reagent set for photometric determination
Range: 0.05–6.00 mg/L Cl₂
Sufficient for: 200 tests
Shelf life: at least 2 years
Sea water suitability: yes
A visual determination is not possible.

VISOCOLOR[®] *ECO* free Chlorine 6 REF 931 219

Type: reagent set for photometric determination
Range: 0.05–6.00 mg/L Cl₂
Sufficient for: 400 tests
Shelf life: at least 2 years
Sea water suitability: yes
A visual determination is not possible.

VISOCOLOR[®] *HE* Chlorine Refill pack

REF 920 015
REF 920 115

Type: highly sensitive test kit
Range: 0.0 · 0.02 · 0.04 · 0.06 · 0.10 · 0.15 · 0.20 · 0.30 · 0.40 · 0.60 mg/L Cl₂
Sufficient for: 160 tests
Shelf life: at least 2 years
Sea water suitability: yes



Test kits for water analysis

Description of individual parameters and tests

Chlorine dioxide

NEW!

ClO₂

Chlorine dioxide is being used as a universal disinfecting reagent, replacing chlorine especially in drinking and swimming pool water processing. Not only does it possess a significantly higher oxidizing power than chlorine but moreover, the formation of toxic THMs (trihalomethanes) is avoided when using chlorine dioxide as the oxidizing agent.

Reaction basis:

Chlorine dioxide reacts, like chlorine, with DPD (*N,N*-diethyl-1,4-phenylenediamine) at a pH value of 5–6 to form a red-violet dye. An additional step makes this test kit even suitable for the determination of chlorine dioxide in presence of chlorine.

VISOCOLOR® ECO Chlorine dioxide

REF 931 021

Refill pack

REF 931 221

Type: colorimetric test kit
Range: <0.2 · 0.2 · 0.4 · 0.6 · 0.8 · 1.1 · 1.7 · 2.3 · 3.8 mg/L ClO₂
Sufficient for: 150 tests
Shelf life: at least 1.5 years
Sea water suitability: no

Chromate

CrO₄²⁻

Chromium compounds can be present in industrial waste waters in trivalent form [chromium(III) ions] and in hexavalent form (chromate and dichromate ions). In order to determine the total chromium all other valencies must be oxidized to chromium(VI). Each test kit is equipped with a detailed description.

Reaction basis:

In sulfuric acid chromate ions react with diphenylcarbazide to form a red-violet dye. (Reaction basis according to DIN EN ISO 7393 G4-2).

VISOCOLOR® ECO Chromium(VI)

REF 931 020

Refill pack

REF 931 220

Type: colorimetric test kit
Range: 0.02 · 0.05 · 0.10 · 0.15 · 0.20 · 0.30 · 0.40 · 0.50 mg/L Cr(VI)
Sufficient for: 140 tests
Shelf life: at least 1.5 years
Sea water suitability: yes



Copper

Cu²⁺

In water, copper(II) can be found in dissolved as well as in undissolved form. Copper(I) compounds and undissolved copper(II) compounds are not determined unless they are decomposed with concentrated nitric acid prior to the test.

Reaction basis:

In weakly alkaline medium copper(II) ions react with cuprizone forming a blue color complex.

VISOCOLOR® ECO Copper

REF 931 037

Refill pack

REF 931 237

Type: colorimetric test kit
Range: 0 · 0.1 · 0.2 · 0.3 · 0.5 · 0.7 · 1.0 · 1.5 mg/L Cu²⁺
Sufficient for: 100 tests
Shelf life: at least 2 years
Sea water suitability: yes

VISOCOLOR® HE Copper

REF 920 050

Refill pack

REF 920 150

Type: highly sensitive test kit
Range: 0.0 · 0.04 · 0.07 · 0.10 · 0.15 · 0.20 · 0.25 · 0.30 · 0.40 · 0.50 mg/L Cu²⁺
Sufficient for: 100 tests
Shelf life: at least 2 years
Sea water suitability: yes



Test kits for water analysis

Description of individual parameters and tests

Cyanide

CN⁻

Cyanide ions are very toxic because they block the iron of the respiratory enzyme and thus inhibit oxygen transport. For humans 1 mg cyanide per kg body weight is considered lethal.

Reaction basis:

Cyanide ions react with chlorine to form cyanogen chloride, which then opens a pyridine ring to form glutacetaldehyde. By aldol condensation with barbituric acid, a violet polymethine dye is produced.

This test detects free cyanide and cyanide complexes which can be destroyed with chlorine. If interfering substances like heavy metal complexes, thiocyanate, sulfide, dyes or aromatic amines are present, a distillation according to DIN 38 405-D 13-2-2 must precede the cyanide test.

For the determination of easily released and total cyanide as well as for the determination of cyanide in stone fruit spirits, please contact MACHEREY-NAGEL for special instructions.

VISOCOLOR® ECO Cyanide Refill pack

REF 931 022
REF 931 222

Type: colorimetric test kit
Range: 0 · 0.01 · 0.02 · 0.03 · 0.05 · 0.07 ·
0.10 · 0.15 · 0.20 mg/L CN⁻
Sufficient for: 100 tests
Shelf life: at least 1 year
Sea water suitability: yes, after dilution (1+3)



VISOCOLOR® HE Cyanide Refill pack

REF 920 028
REF 920 128

Type: highly sensitive test kit
Range: 0.0 · 0.002 · 0.004 · 0.007 · 0.010 ·
0.015 · 0.020 · 0.025 · 0.030 ·
0.040 mg/L CN⁻
Sufficient for: 55 tests
Shelf life: at least 1 year
Sea water suitability: yes

Cyanuric acid

Cya

Chlorine for disinfection in swimming pools is degraded by intensive UV radiation. A common stabilizer for chlorine in swimming pools is cyanuric acid. Chloroisocyanic acid is also used directly as a disinfecting agent. Due to legal regulations, as good and careful pool managers, many private and public owners use rapid tests for frequent monitoring of cyanuric acid.

Reaction basis:

Turbidity measurement

Cyanuric acid forms a fine precipitate with a triazine derivative.

The turbidity caused by this reaction can be measured visually or photometrically to determine the cyanuric acid concentration.

Turbidities interfere and must be filtered prior to the analysis.

VISOCOLOR® ECO Cyanuric acid Refill pack

REF 931 023
REF 931 223

Type: colorimetric test kit
Range: 10 · 15 · 20 · 30 · 40 · 60 · 80 ·
100 mg/L Cya
Sufficient for: 100 tests
Shelf life: at least 1.5 years
Sea water suitability: yes

DEHA (Diethylhydroxylamine)

DEHA

In boiler feed water treatments, the carcinogenic hydrazine is more and more replaced by diethylhydroxylamine (DEHA) to remove oxygen.

Reaction basis:

Measurement of the reduction properties of DEHA for iron(III) ions and determination of the iron(II) ions formed. Strictly observe the temperature and reaction time since they strongly influence the color intensity.

Iron(II) ions interfere. Modified test procedure eliminate this source of error.

VISOCOLOR® ECO DEHA

REF 931 024

Refill pack

REF 931 224

Type: colorimetric test kit
Range: 0 · 0.01 · 0.03 · 0.05 · 0.10 · 0.15 ·
0.20 · 0.25 · 0.30 mg/L DEHA
Sufficient for: 125 tests
Shelf life: at least 1 year
Sea water suitability: yes



Test kits for water analysis

Description of individual parameters and tests

Dithionite



Dithionites ($\text{S}_2\text{O}_4^{2-}$), above all sodium dithionite (so-called hydrosulfite), are important aids in dyeing, textile and paper industries due to their reducing properties. They are especially used for vat dyeing, for bleaching of wood-based paper, sugar, sirup, gelatine, starche, molasses, saccharine juice, soap, technical fats, as decolorant for textiles, for desilvering of used fixing baths etc. Dithionite can be determined using a special method with *VISOCOLOR*[®] Sulfite SU 100. Please contact MACHEREY-NAGEL for special instructions.

Fluoride



Normally, the content of fluoride in surface and ground water is lower than 1 mg/L. Permanent consumption of water containing more than 2 mg/L fluoride can cause a drinking water fluorese (a stained discoloration of the tooth enamel). On the other hand a content of fluoride lower than 0.5 mg/L could lead to an increased cavity risk. The optimal content of fluoride in drinking water is about 1 mg/L. To avoid cavity problems, in some countries the content of fluoride in drinking water is increased artificially. However, WHO and EC drinking water regulations recommend a threshold value of 1.5 mg/L F^- .

In some European countries, the threshold value for mineral water is also 1.5 mg/L F^- . For water used in the preparation of baby food a threshold value of 0.7 mg/L F^- is recommended. Natural mineral waters with a fluoride content of more than 5 mg/L must be labelled with a warning.

Reaction basis:

Photometric determination of fluoride with 1,8-dihydroxy-2-(4-sulfophenylazo)naphthalene-3,6-disulfonic acid (SPADNS) with the photometer PF-12.

VISOCOLOR[®] ECO Fluoride

REF 931 227

Type: reagent set for photometric determination
Range: 0.1–2.0 mg/L F^-
Sufficient for: 150 tests
Shelf life: at least 1.5 years
Sea water suitability: yes, after distillation

A visual determination is not possible

Hardness (total and residual)



The total hardness of water is based on its content of alkaline earth ions (calcium and magnesium ions). This content depends on the geological conditions and may vary considerably. Knowledge of the total hardness is important for the use of water in industrial as well as municipal applications, e.g. in the household as wash water or as boiler feed water in industry.

Reaction bases:

(a) Complexometric titration
in accordance with DIN 38406 E3 and DIN 38409 H6.

(b) Colorimetry with a mixed indicator
Copper(II) ions can delay or (in higher concentrations) even block the color change of the indicator. For this reason allow enough water to run through copper pipes prior to sampling.

VISOCOLOR[®] alpha total Hardness

REF 935 042

Type: titration test kit
Range: 1 drop \triangleq 1.25 $^{\circ}\text{e} \triangleq$ 17.8 mg/L CaCO_3
Reaction basis: (a) titration
Sufficient for: 100 tests with an average hardness of 12.5 $^{\circ}\text{e}$
Shelf life: at least 1.5 years
Sea water suitability: yes, after dilution (1:30)

VISOCOLOR[®] ECO total Hardness

REF 931 029

Type: titration test kit
Range: 1 drop \triangleq 1.25 $^{\circ}\text{e} \triangleq$ 17.8 mg/L CaCO_3
Reaction basis: (a) titration
Sufficient for: 110 tests with an average hardness of 12.5 $^{\circ}\text{e}$
Shelf life: at least 1.5 years
Sea water suitability: yes, after dilution (1:30)



VISOCOLOR[®] HE total Hardness H 20 F Refill pack

REF 915 005
REF 915 205

Type: titration test kit
Range: 0.6–25.0 $^{\circ}\text{e}$ or
0.1–3.6 mmol/L Ca^{2+}
1 gradation mark = 0.5 $^{\circ}\text{d}$ =
0.1 mmol/L Ca^{2+}
Reaction basis: (a) titration
Sufficient for: about 200 tests with an average hardness of 12.5 $^{\circ}\text{e}$ or 1.8 mmol/L Ca^{2+}
Shelf life: at least 1.5 years
Sea water suitability: yes, after dilution (1:30)

Test kits for water analysis

Description of individual parameters and tests

VISOCOLOR® alpha residual Hardness REF 935 080

Type: colorimetric test kit
 Range: 0.00 · 0.05 · 0.10 · 0.19 · 0.38 °e
 Reaction basis: (b) colorimetry
 Sufficient for: 200 tests
 Shelf life: at least 1 year
 Sea water suitability: no

VISOCOLOR® HE Hardness H 2 REF 915 002 **Refill pack** REF 915 202

Type: titration test kit
 Range: 0.06–2.50 °e or
 0.01–0.36 mmol/L Ca²⁺
 1 gradation mark = 0.06 °e or
 0.01 mmol/L Ca²⁺
 Reaction basis: (a) titration
 Sufficient for: 200 tests with an average hardness
 of 1.25 °e or 0.18 mmol/L Ca²⁺
 Shelf life: at least 1.5 years
 Sea water suitability: no

VISOCOLOR® ECO additive reagent Z-1 REF 931 929

to eliminate copper ions during determination of total hardness

Hydrazine



Hydrazine is used to destroy residual oxygen in boiler feed water and condensate water, for example in power plants, to avoid corrosion of the boiler casing. Reaction products are merely nitrogen gas and water, thus the salt load of the water is kept low.

Because of its highly reactive properties, hydrazine is also used as fuel in aviation and astronautics.

Hydrazine is toxic and has a highly toxic effect on water organisms. Hydrazine can be absorbed through the skin. Therefore, water and waste water with potential content of hydrazine must be monitored and tested.

Reaction basis:

DIN method: In acidic solution hydrazine reacts with 4-dimethylaminobenzaldehyde to form a yellow/orange colored compound (Reaction basis according to DIN 38413-P1).

VISOCOLOR® ECO Hydrazine REF 931 030 **Refill pack** REF 931 230

Type: titration test kit
 Range: 0 · 0.05 · 0.10 · 0.15 · 0.20 · 0.25 ·
 0.30 · 0.40 mg/L N₂H₄
 Sufficient for: 130 tests
 Shelf life: at least 1 year
 Sea water suitability: yes

Hydrosulfite



see Dithionite, page 62

Iron

Fe

Natural waters as well as waste waters often contain some iron. It can be present as Fe(II) or Fe(III) ions if the pH value of the water is below 3 or if the water is free of oxygen. At higher pH values Fe(III) forms an insoluble oxyhydrate. Frequently the sample is in a transition state with finely dispersed iron oxides. In waste water and natural water containing humic acids iron is often present in the form of a complex salt. The VISOCOLOR® test kits only determine iron which is present in dissolved form as Fe²⁺ or Fe³⁺ ions. Iron complexes are not covered unless they are decomposed by oxidation with nitric acid and sulfuric acid.

Reaction basis:

Triazin-Methode: iron(II) ions react with a triazine derivative to form a violet complex. Iron(III) ions are reduced and thus also determined.

VISOCOLOR® ECO Iron 1 **Refill pack**

NEW!

REF 931 025
 REF 931 225

Type: colorimetric test kit
 Range: 0 · 0.04 · 0.07 · 0.10 · 0.15 · 0.20 ·
 0.30 · 0.50 · 1.0 mg/L Fe
 Sufficient for: 200 tests
 Shelf life: at least 2 years
 Sea water suitability: yes

VISOCOLOR® ECO Iron 2 **Refill pack**

REF 931 026
 REF 931 226

Type: colorimetric test kit
 Range: 0 · 0.04 · 0.07 · 0.10 · 0.15 · 0.20 ·
 0.30 · 0.50 · 1.0 mg/L Fe
 Sufficient for: 100 tests
 Shelf life: at least 2 years
 Sea water suitability: yes

VISOCOLOR® HE Iron **Refill pack**

REF 920 040
 REF 920 140

Type: highly sensitive test kit
 Range: 0.0 · 0.01 · 0.02 · 0.03 · 0.04 · 0.05 ·
 0.07 · 0.10 · 0.15 · 0.20 mg/L Fe
 Sufficient for: 300 tests
 Shelf life: at least 2 years
 Sea water suitability: no



Test kits for water analysis

Description of individual parameters and tests

Magnesium

Mg²⁺

In order to differentiate between the hardness constituents calcium and magnesium one can determine the total hardness (see hardness, page 62) and the calcium hardness (see calcium, page 57). The difference of both is the magnesium content – an important parameter in food and building industries.

Manganese

Mn

Natural waters contain manganese in the divalent, soluble form as well as in colloidal tri- and tetravalent states. Valencies are interchangeable due to oxidation-reduction reactions taking place in the water. The test procedure determines all oxidation states of manganese.

Reaction basis:

In alkaline solution manganese ions react with formaldoxime to form an orange-red complex. (Reaction basis according to DIN 38406-E2).

VISOCOLOR® ECO Manganese

REF 931 038

Refill pack

REF 931 238

Type: colorimetric test kit
Range: 0 · 0.1 · 0.2 · 0.3 · 0.5 · 0.7 · 0.9 · 1.2 · 1.5 mg/L Mn

Sufficient for: 70 tests
Shelf life: at least 1.5 years
Sea water suitability: yes

VISOCOLOR® HE Manganese

REF 920 055

Refill pack

REF 920 155

Type: highly sensitive test kit
Range: 0.0 · 0.03 · 0.06 · 0.10 · 0.15 · 0.20 · 0.25 · 0.30 · 0.40 · 0.50 mg/L Mn

Sufficient for: 100 tests
Shelf life: at least 1.5 years
Sea water suitability: no



Nickel

Ni²⁺

Nickel can be present in industrial waste water. It occurs as the divalent ion or as a nickel complex.

Reaction basis:

In ammonia solution nickel ions react with diacetyldioxime after oxidation with bromine to form a reddish-brown dye. Insoluble nickel compounds (e.g. nickel cyanide, nickel carbonate) and the nickel cyano complexes are not determined.

VISOCOLOR® ECO Nickel

REF 931 040

Refill pack

REF 931 240

Type: colorimetric test kit
Range: 0 · 0.1 · 0.2 · 0.3 · 0.5 · 0.7 · 0.9 · 1.2 · 1.5 mg/L Ni²⁺

Sufficient for: 150 tests
Shelf life: at least 1.5 years
Sea water suitability: yes, after dilution (1+9)



Nitrate

NO₃⁻

Nitrates occur in most ground and surface waters in concentrations of up to 20 mg/L. In addition to mere geological influences, nitrate concentrations can also increase due to agricultural sources (fertilisers). The VISOCOLOR® test kits are intended for the determination of nitrate in surface and drinking water and in industrial waste waters which do not contain high concentrations of interfering ions.

Reaction basis:

Nitrate is reduced to nitrite with an inorganic reducing agent. Nitrite is then diazotised with an aromatic amine and simultaneously coupled to form an azo dye.

VISOCOLOR® alpha Nitrate

REF 935 065

Type: colorimetric test kit
Range: 2 · 8 · 15 · 30 · 50 mg/L NO₃⁻
Sufficient for: 100 tests
Shelf life: at least 1.5 years
Sea water suitability: yes

VISOCOLOR® ECO Nitrate

REF 931 041

Refill pack

REF 931 241

Type: colorimetric test kit
Range: 0 · 1 · 3 · 5 · 10 · 20 · 30 · 50 · 70 · 90 · 120 mg/L NO₃⁻

Sufficient for: 110 tests
Shelf life: at least 1.5 years
Sea water suitability: yes

Test kits for water analysis

Description of individual parameters and tests

Nitrite



In surface waters nitrite ions are generally present in low concentrations. Their presence in ground water is less common. In waste waters nitrite frequently occurs, even in fairly high concentrations.

Reaction basis:

Sulfanilamide/sulfanilic acid is diazotized by nitrite in acidic solution. The diazonium salt is coupled with a naphthylamine to form an intensively colored azo dye.

VISOCOLOR® alpha Nitrite

REF 935 066

Type: colorimetric test kit
Range: 0.05 · 0.10 · 0.25 · 0.5 · 1.0 mg/L NO_2^-
Sufficient for: 200 tests
Shelf life: at least 1.5 years
Sea water suitability: yes

VISOCOLOR® ECO Nitrite

REF 931 044

Refill pack

REF 931 244

Type: colorimetric test kit
Range: 0 · 0.02 · 0.03 · 0.05 · 0.07 · 0.1 · 0.2 · 0.3 · 0.5 mg/L NO_2^-
Sufficient for: 120 tests
Shelf life: at least 1.5 years
Sea water suitability: yes



VISOCOLOR® HE Nitrite

REF 920 063

Refill pack

REF 920 163

Type: highly sensitive test kit
Range: 0.0 · 0.005 · 0.010 · 0.015 · 0.02 · 0.03 · 0.04 · 0.06 · 0.08 · 0.10 mg/L NO_2^-
Sufficient for: 150 tests
Shelf life: at least 2 years
Sea water suitability: yes

Oxygen



The solubility of oxygen in water depends on the temperature, the pressure and other water components. The oxygen content of water at the time of sampling is often quoted as the percentage of possible saturation.

Reaction basis:

Oxygen determination according to Winkler: in alkaline solution dissolved oxygen oxidises manganese(II) ions to higher-valent manganese hydroxides. In a strongly acidic medium these form manganese(III) ions which can be determined either titrimetrically or colorimetrically after addition of a color reagent.

VISOCOLOR® ECO Oxygen

REF 931 088

Refill pack

REF 931 288

Type: colorimetric test kit
Range: 0 · 1 · 2 · 3 · 4 · 6 · 8 · 10 mg/L O_2
Sufficient for: 50 tests
Shelf life: at least 1 year
Sea water suitability: yes

When ordering this test kit for the first time, you also need an oxygen bottle, REF 915 498

VISOCOLOR® HE Oxygen SA 10

REF 915 009

Refill pack

REF 915 209

Type: titration test kit, according to DIN EN 25 813
Range: 0.2–10 mg/L O_2
1 gradation mark = 0.2 mg/L
Sufficient for: about 100 tests with an average oxygen content of 9 mg/L
Shelf life: at least 1.5 years
Sea water suitability: yes



In combination with the BOD₅ accessories package (REF 916 918) and the BOD₅ nutrient mixture without N-allylthiourea (ATU) (REF 918 994) or the BOD₅ nutrient mixture Plus with ATU (REF 918 995) this test kit can also be used for determination of the BOD₅ (sufficient for 25–50 samples). Preparation of the samples is performed using the so-called dilution principle according to DIN ISO 1899-1-H51.

Test kits for water analysis

Description of individual parameters and tests

pH value

The pH value indicates whether a water reacts acidic, alkaline or neutral. It is determined by the concentration of hydrogen ions. All biological processes in water are tied to specific pH ranges. For municipal and industrial applications, too, the control of specific pH limits is important since e.g. the efficiency of sewage plants or the corrosive potential of tap water on pipes depend in the pH value. Contrary to pH indicator papers, **VISOCOLOR®** pH test kits can also be used to determine accurate pH values in unbuffered sample solutions.

Reaction basis:

A special mixture of indicator dyes shows a characteristic color for every pH value within the range of the kit.

The favorable ratio between sample volume and amount of indicator minimizes the indicator error (acid-base-error). This allows reliable pH measurements in weakly buffered solutions as well. High contents of neutral salts and colloids as well as organic solvent concentrations above 10 % can cause wrong results.

VISOCOLOR® alpha pH 5–9 REF 935 075

Type: colorimetric test kit
Range: pH 5.0 · 5.5 · 6.0 · 6.5 · 7.0 · 7.5 · 8.0 · 8.5 · 9.0
Sufficient for: 200 tests
Shelf life: at least 3 years
Sea water suitability: yes

VISOCOLOR® ECO pH 4.0–9.0 REF 931 066 Refill pack REF 931 266

Type: colorimetric test kit
Range: pH 4.0 · 5.0 · 6.0 · 6.5 · 7.0 · 7.5 · 8.0 · 8.5 · 9.0
Sufficient for: 450 tests
Shelf life: at least 3 years
Sea water suitability: yes

VISOCOLOR® ECO pH 6.0–8.2 REF 931 270

Type: reagent set for photometric determination
Range: pH 6.0–8.2
Sufficient for: 150 tests
Shelf life: at least 1.5 years
Sea water suitability: yes

A visual colorimetric determination is not possible.

VISOCOLOR® HE pH 4.0–10.0 REF 920 074 Refill pack REF 920 174

Type: highly sensitive test kit
Range: pH 4.0 · 5.0 · 5.5 · 6.0 · 6.5 · 7.0 · 7.5 · 8.0 · 8.5 · 9.0 · 10.0
Sufficient for: 500 tests
Shelf life: at least 2 years
Sea water suitability: yes

For the determination of the pH value, see also Swimming pool page 69

pH

Phosphate

PO₄³⁻

The phosphate content in surface water has direct consequences for its ability to support the growth of certain organisms. Since increasing quantities of phosphates are being fed through domestic waste waters into rivers and lakes, these waters have a tendency towards eutrophication. Exact data about the phosphate content are important for boiler water and feed water. Precise dosing of phosphates in these waters can inhibit formation of boiler scale. Pyro-, meta- and poly-phosphates are not determined with **VISOCOLOR®** test kits. The determination of total phosphate requires a decomposition prior to the test.

Reaction basis:

Ammonium molybdate reacts with phosphate ions to form phosphomolybdic acid which is reduced to molybdenum blue (Reaction basis according to DIN EN ISO 6878-D11).

VISOCOLOR® alpha Phosphate REF 935 079

Type: colorimetric test kit
Range: 2 · 5 · 10 · 15 · 20 mg/L PO₄³⁻
Sufficient for: 70 tests
Shelf life: at least 2 years
Sea water suitability: yes

VISOCOLOR® ECO Phosphate REF 931 084 Refill pack REF 931 284

Type: colorimetric test kit
Range: 0 · 0.2 · 0.3 · 0.5 · 0.7 · 1 · 2 · 3 · 5 mg/L PO₄-P
Sufficient for: 90 tests
Shelf life: at least 3 years
Sea water suitability: yes



VISOCOLOR® HE Phosphate REF 920 082 Refill pack REF 920 182

Type: highly sensitive test kit
Range: 0.0 · 0.05 · 0.10 · 0.15 · 0.20 · 0.3 · 0.4 · 0.6 · 0.8 · 1.0 mg/L P
Sufficient for: 300 tests
Shelf life: at least 2 years
Sea water suitability: yes

Test kits for water analysis

Description of individual parameters and tests

VISOCOLOR® HE Phosphate (DEV)

REF 920 080

Refill pack

REF 920 180

Type: highly sensitive test kit
 Range: 0.0 · 0.01 · 0.02 · 0.03 · 0.05 · 0.07 · 0.10 · 0.15 · 0.20 · 0.25 mg/L P
 Sufficient for: 100 tests
 Shelf life: at least 2 years
 Sea water suitability: yes



Phosphonate

PO_x

Phosphonates are used as complexing agents for softening process and cooling water. They can be determined with NANOCOLOR® NanOx Metal and VISOCOLOR® ECO Phosphate. Please ask for our special instructions.

Potassium

K⁺

The natural potassium content in ground water is generally about 1–2 mg/L K⁺. Higher values may indicate faecal contaminations, but can also originate from potassium fertilizers. For the growth of plants and animals, potassium is an essential factor. Especially in agriculture the determination of potassium therefore gains increasing importance.

Reaction basis:

Potassium reacts with sodium tetraphenylborate to form a precipitate. Under defined conditions this turbidity can be used for concentration measurements.

Turbidities interfere and have to be filtered prior to the test. Good reproducibility is obtained for drinking, surface and ground water. For polluted waste waters low potassium values may be measured.

VISOCOLOR® ECO Potassium

REF 931 032

Refill pack

REF 931 232

Type: turbidity test kit
 Range: 2 · 3 · 4 · 6 · 8 · 10 · 15 mg/L K⁺
 Sufficient for: 60 tests
 Shelf life: at least 3 years
 Sea water suitability: yes, after dilution (1+1)

Residual hardness

°e

see hardness (total and residual), page 62

Silica/silicon

SiO₂

Natural water contains silica in different amounts, depending on the geological conditions. The silica occurs partly as soluble silicate, partly colloidal as polysilicic acids. The silica content in water for medium and high pressure boilers may not exceed certain limits. Therefore, in power plants constant measurement of the silica content is required.

Reaction basis:

In acidic solution soluble silica or silicates react with ammonium molybdate to form yellow silicomolybdic acid, which is reduced to silico-molybdenum blue with a reducing agent. (Reaction basis according to DIN EN ISO 16264-H57).

VISOCOLOR® ECO Silica

REF 931 033

Refill pack

REF 931 233

Type: colorimetric test kit
 Range: 0 · 0.2 · 0.4 · 0.6 · 1.0 · 1.5 · 2.0 · 2.5 · 3.0 mg/L SiO₂

Sufficient for: 80 tests
 Shelf life: at least 3 years
 Sea water suitability: yes

VISOCOLOR® HE Silicon

REF 920 087

Refill pack

REF 920 187

Type: highly sensitive test kit
 Range: 0.0 · 0.01 · 0.02 · 0.03 · 0.05 · 0.07 · 0.10 · 0.15 · 0.20 · 0.30 mg/L Si

Sufficient for: 120 tests
 Shelf life: at least 2 years
 Sea water suitability: yes



Test kits for water analysis

Description of individual parameters and tests

Sodium

 Na^+

About 2.43 % of the uppermost crust of the earth, with a thickness of 16 km, consist of chemically bonded sodium. It is the 6th most abundant element on our planet. In rocks and minerals it is mostly present as silicate (e.g. soda feldspar or albite) and as sodium chloride (rock salt), but also as carbonate, nitrate (especially Chile saltpeter) and cryolite as well as in many other minerals. Each ton of sea water contains an average of 27 kg common salt (10.6 kg sodium), so, sodium makes up for 77 % of all salts present in sea water. Under the assumption, that in natural waters beside sodium mainly calcium and magnesium can occur as cations, it is easily possible to determine the sodium content using *VISOCOLOR*® titration test kits. Please ask for the respective special instructions.

Sulfate

 SO_4^{2-}

The determination of sulfate ions is of special importance to evaluate the aggressiveness of a water towards concrete. Sulfate, a component of natural water, is formed in many manufacturing processes and is found in waste waters even after neutralization.

Reaction basis:

Sulfate reacts with barium ions forming a precipitate of barium sulfate. Under defined conditions this turbidity can be used for concentration measurements.

Turbidities of the sample interfere and have to be filtered. Good reproducibility is obtained for drinking, surface and ground water. For polluted waste waters low concentrations are measured.

VISOCOLOR® ECO Sulfate

Refill pack

Type: turbidity test kit
 Range: 25 · 30 · 35 · 40 · 50 · 60 · 70 · 80 · 100 · 120 · 150 · 200 mg/L SO_4^{2-}
 Sufficient for: 100 tests
 Shelf life: at least 3 years
 Sea water suitability: yes, after dilution (1:50)

REF 931 092

REF 931 292

Sulfide

 S^{2-}

In water sulfides can be present as dissolved hydrogen sulfide or as hydrosulfide or sulfide ions.

Reaction basis:

N,N-Dimethyl-1,4-phenylene diamine reacts with hydrogen sulfide to form an unstable compound which rearranges to leucomethylene blue. Oxidation with iron(III) ions yields methylene blue.

VISOCOLOR® ECO Sulfide

Refill pack

Type: colorimetric test kit
 Range: 0.1 · 0.2 · 0.3 · 0.4 · 0.5 · 0.6 · 0.7 · 0.8 mg/L S^{2-}
 Sufficient for: 90 tests
 Shelf life: at least 3 years
 Sea water suitability: yes

REF 931 094

REF 931 294



Sulfite

 SO_3^{2-}

Sulfite ions are not present in natural, unpolluted waters, however they often appear in large quantities in industrial waste waters (e.g. paper mills, dye works).

Reaction basis:

Titrimetric determination by addition of iodine solution and reverse titration of unused iodine with sodium thiosulfate.

Oxidizing and reducing substances interfere.

VISOCOLOR® ECO Sulfite

Type: titration test kit
 Range: 1 drop = 1 mg/L SO_3^{2-}
 Sufficient for: 60 tests with an average Sulfite concentration of 10 mg/L SO_3^{2-}
 Shelf life: at least 1 year
 Sea water suitability: yes

REF 931 095

VISOCOLOR® HE Sulfite SU 100

Refill pack

Type: titration test kit
 Range: 2–100 mg/L SO_3^{2-}
 1 gradation mark = 2 mg/L
 Sufficient for: 100 tests with an average Sulfite content of 100 mg/L SO_3^{2-}
 Shelf life: at least 3 years
 Sea water suitability: yes

REF 915 008

REF 915 208

Test kits for water analysis

Description of individual parameters and tests

Swimming pool

Cl₂ pH

Chlorination of water in swimming pools with chlorine or chlorine compounds takes place within certain limits. One distinguishes between free chlorine and bound chlorine (chloramines); the sum of both is called total chlorine. The content of free chlorine should be between 0.3 and 0.6 mg/L. Chlorination alters the pH value of the swimming pool water. The ideal pH value is pH 7.4. This prevents formation of malodorous pollutants and irritants to the mucous membrane, and prevents damage to the water treatment system.

Reaction basis:

Free chlorine reacts with *N,N*-diethyl-1,4-phenylene diamine (DPD) to form a red-violet dye. The pH value is determined with phenol red as indicator.

For test kits for the determination of chlorine alone see VISOCOLOR® Chlorine (see page 59)

VISOCOLOR® ECO Swimming pool Refill pack

REF 931 090
REF 931 290

Type: colorimetric test kit
Range: < 0.1 · 0.2 · 0.3 · 0.4 · 0.6 · 0.9 · 1.2 · 2.0 mg/L Cl₂
pH 6.9 · 7.2 · 7.4 · 7.6 · 7.8 · 8.2
Sufficient for: 150 tests each
Shelf life: at least 1.5 years
Sea water suitability: yes

Zinc

Zn²⁺

Zinc is one of the most often used metals for surface finishing. Its content in waste water from e.g. electroplating companies has to be monitored regularly.

Reaction basis:

At a pH of 8.5–9.5 zinc ions react with zincon to form a color complex. Acidic, alkaline and buffered samples have to be adjusted to pH 9 for the test.

VISOCOLOR® ECO Zinc Refill pack

REF 931 098
REF 931 298

Type: colorimetric test kit
Range: 0 · 0.5 · 1 · 2 · 3 mg/L Zn²⁺
Sufficient for: 120 tests
Shelf life: at least 1 year
Sea water suitability: yes, after dilution (1+9)



visocolor®

Reagent cases

Combinations of different **VISOCOLOR®** test kits

Reagent cases for all purposes

- Rugged cases with premium foam inlays
- Pre-packed or empty cases for individual solutions
- With and without photometer PF-12

MACHEREY-NAGEL reagent cases are flexible tools for all areas of water and soil analysis. The cases are ideal for testing water bodies, for fish farming and a wide variety of other

applications within the area of environmental management and water monitoring.

The reagent cases enable users to get results fast and directly at the point of interest. Previous chemical knowledge or experience is not required to run any of the tests or to use the cases effectively.

The cases include all necessary test instructions and analytical accessories for especially easy and convenient handling.

Ordering information

Type	For the determination of	REF
Reagent case for soil analysis		
VISOCOLOR® reagent case for soil analysis*	Phosphate, potassium, nitrogen (ammonium, nitrite and nitrate), pH, soil structure, incl. all requisite equipment and accessories	931 601
Reagent case for schools		
VISOCOLOR® SCHOOL *	Ammonium, nitrate, nitrite, phosphate, pH, total hardness	933 100
Reagent case for water analysis · without photometer		
VISOCOLOR® ECO reagent case*	Ammonium, carbonate hardness, total hardness, nitrate, nitrite, pH, phosphate	931 301
VISOCOLOR® ECO reagent case, empty	For individual combination of up to 8 test kits of type VISOCOLOR® ECO	931 303
VISOCOLOR® reagent case*	Alkalinity, ammonium, total hardness, nitrite, pH, phosphate, oxygen, temperature	931 304
VISOCOLOR® reagent case, empty	For combination with products from the Rapid Tests and the VISOCOLOR® lines	931 305
Reagent case for water analysis · with photometer PF-12		
Reagent case "Environmental analysis"*	Ammonium, carbonate hardness, iron, total hardness, nitrate, nitrite, pH, phosphate, incl. photometer PF-12	914 303
VISOCOLOR® reagent case with PF-12	For combination with products from the Rapid Tests and the VISOCOLOR® lines	914 301
* This product contains harmful substances which must be specially labeled as hazardous. For detailed information please see MSDS.		

Soil analysis made easy

- Easy sample preparation
- Fast evaluation of nutrients in soil
- Convenient and reliable determination even in the field

VISOCOLOR® reagent case for soil analysis REF 931 601

Thorough analysis is a corner stone to support and maintain healthy, productive and biologically active soil. To effectively and efficiently plan all measures that affect the soil (fertilization, liming, etc.) it is crucial to determine the important soil parameters first. The **VISOCOLOR®** reagent case for soil analysis is the perfect companion for economical, fast and convenient soil analysis, both in the field or in your laboratory, including the following tests:

QUANTOFIX® test strips:

- Ammonium 10–400 mg/L NH_4^+
- Nitrate / Nitrite 10–500 mg/L NO_3^- / 1–80 mg/L NO_2^-

VISOCOLOR® ECO test kits:

- Potassium 2–15 mg/L K^+

VISOCOLOR® HE test kits:

- Phosphate 1–20 mg P/100 g
- pH 4.0–10.0

pH-Fix test strip

- pH 2.0–9.0

New design!



Depending on the parameter to be determined, the soil extracts are prepared with Calcium-Acetate-Lactate (CAL) or with CaCl_2 solution. If required by national regulations or because of local geological conditions, the **VISOCOLOR®** reagent case for soil analysis can also be used with other extraction solutions.

Soil analysis in no time


- Manual with background information
- Ideal combination of test papers and colorimetric test kits
- Incl. additional analytical tools, such as scale, sieve etc.

In addition to the reagent case, various **VISOCOLOR®** reagents can also be used for soil analysis. Please request information about suitable extraction solutions.

Combinations of different VISOCOLOR® test kits

Water analysis – it's that simple

- Easy test procedures
- Safe results
- Flexible applications
- Incredible value



NEW!



VISOCOLOR® SCHOOL reagent case REF 933 100

The VISOCOLOR® SCHOOL reagent case is especially designed for schools and caters to the needs of both students and teachers. All reagents are approved to be used in schools in Germany (GUV-SR 2004 directive) and can be disposed of easily just down the drain without any harm to the environment.

VISOCOLOR® SCHOOL reagents:


Ammonium	0.2–3 mg/L NH_4^+
Nitrate	1–90 mg/L NO_3^-
Nitrite	0.02–0.5 mg/L NO_2^-
Phosphate	0.5–15 mg/L PO_4^{3-}
pH	4.0–9.0
Total hardness	1 drop \triangleq 17.8 mg/L CaCO_3

The newly developed tests are based on our proven VISOCOLOR® system and were developed as high quality VISOCOLOR® SCHOOL test specifically for schools.

Great features – incredible price

- Maximum safety due to exact labeling of all reagents
- Safe results using color and turbidity compensation
- High sensitivity down to the values of German drinking water standards
- Safe for the environment and easy disposal of used tests

High accuracy




- Color matching with all-in-one color scale
- Compensation of turbidity and colors
- Reaction principles based on international standards

Ordering information

Type	REF
VISOCOLOR® SCHOOL refill pack*	933 200
VISOCOLOR® SCHOOL color chart	933 300
VISOCOLOR® SCHOOL manual (German / English)	933 150
* This product contains harmful substances which must be specially labeled as hazardous. For detailed information please see MSDS.	

Simple and safe



- Easy test procedures with pictogram instructions
- Color-coded reagent bottles
- No additional accessories required

Extensive manual



- Background information about the parameters
- Declaration of the reaction basis as well as reaction equation
- Data about interferences

Optionally, the case can be equipped with a pack of pH-Fix or QUANTOFIX® test strips as well.

Reagent cases

Combinations of different **VISOCOLOR®** test kits

Reagent cases for water analysis without photometer

- Easy to handle mini-lab
- For exact and accurate determination
- Expandable with one additional **VISOCOLOR® ECO** test
- Compact case for extensive analysis
- Highly sensitive **VISOCOLOR® HE** tests included
- Infinite options

VISOCOLOR® ECO reagent case

REF 931 301

For mobile use with test instructions for visual determination, analytical accessories and the following tests:

VISOCOLOR® ECO test kits:

Ammonium	0.2–3 mg/L NH_4^+
Carbonate Hardness	1 drop \triangleq 1.25 °e
total Hardness	1 drop \triangleq 1.25 °e
Nitrate	1–120 mg/L NO_3^-
Nitrite	0.02–0.5 mg/L NO_2^-
pH	pH 4.0–9.0
Phosphate	0.2–5 mg/L P



VISOCOLOR® reagent case

REF 931 304

For mobile use with test instructions for colorimetric and titrimetric determinations with **VISOCOLOR®** tests, analytical accessories as well as the following test kits:

VISOCOLOR® ECO test kits:

Ammonium 3	0.2–3 mg/L NH_4^+
Nitrite	0.02–0.5 mg/L NO_2^-
Phosphate	0.2–5 mg/L P
pH 4.0–9.0	pH 4.0–9.0

VISOCOLOR® HE test kits:

Alkalinity AL 7	0.2–7.0 mmol/L
total Hardness H 20 F	0.6–25.0 °e
Oxygen SA 10	0.2–10.0 mg/L O_2



VISOCOLOR® ECO reagent case, empty

REF 931 303

The empty **VISOCOLOR® ECO** reagent case allows the case to be equipped in an individual and highly flexible way with 8 different **VISOCOLOR® ECO** test kits.

In addition to the colorimetric and titrimetric **VISOCOLOR® ECO** tests the reagent case can be equipped with **VISOCOLOR® ECO** turbidity tests for potassium, cyanuric acid and sulphate as well. Furthermore the cases can be equipped with the **VISOCOLOR® ECO** Oxygen test including the necessary oxygen flask.

VISOCOLOR® reagent case, empty

REF 931 305

The empty **VISOCOLOR®** reagent case enable you to freely combine every combination of **VISOCOLOR®** tests, pH indicator papers, pH-Fix test strips, qualitative test papers and semi-quantitative **QUANTOFIX®** test strips.

Additionally, the case can be equipped with **VISOCOLOR® ECO** turbidity tests for potassium, cyanuric acid and sulfate as well as with the **VISOCOLOR® ECO** or **VISOCOLOR® HE** Oxygen test including oxygen flask.

Reagent cases

Combinations of different *VISOCOLOR*® test kits

Reagent cases with photometer PF-12

- The all-rounder for any purposes
- Increased accuracy and reproducibility
- Infinite options

VISOCOLOR® reagent case “Environmental analysis”

REF 914 303

Reagent case with photometer PF-12 (see page 74) incl. instrument manual, *VISOCOLOR*® ECO test instructions as well as test instructions for titration test kits, analytical accessories and the following test kits:

VISOCOLOR® ECO test kits:

Ammonium 15	0.5–15 mg/L NH_4^+
Iron	0.04–1.0 mg/L Fe
Nitrate	1–120 mg/L NO_3^-
Nitrite	0.02–0.5 mg/L NO_2^-
Phosphate	0.2–5 mg/L P
pH 4.0–9.0	pH 4.0–9.0

VISOCOLOR® HE test kits:

Carbonate hardness C 20	0.6–25.0 °e
total Hardness H 20 F	0.6–25.0 °e



Customized case solutions

Catering to individual customer needs is of the great importance to MACHEREY-NAGEL. Even though our case solutions already provide a high level of flexibility, we recognize that some customers may have specific requirements outside our existing case solutions. Therefore, we offer entirely individual solutions with a foam inlay designed exactly to your specifications and testing needs. Starting at a minimum of 50 cases, we can provide you with a case that perfectly fits your personal requirements. In addition to a special inlay, we also offer readily packed cases starting at a minimum quantity of 50 cases as well.

Thus, within our highly flexible case range, we can provide virtually any customer with the perfect testing and transportation solution.

For a case solution for test strip and test papers without *VISOCOLOR*® test kits, see page 46.



VISOCOLOR® reagent case with PF-12, without test kits

REF 914 301

The premium foam inlays of the empty *VISOCOLOR*® reagent case with PF-12 allow the case to be equipped in an individual and highly flexible way with every combination of the *VISOCOLOR*® tests as well as pH indicator papers, pH-Fix test strips, qualitative test papers and semi-quantitative QUANTOFIX® test strips.

visocolor®

Photometer for water analysis

Compact photometer PF-12

Maximum flexibility


- More than 100 preprogrammed methods
- Automatic wavelength adjustment
- Backlit graphic display with intuitive userguidance
- Data storage according to GLP
- USB interface for data transfer, update and power supply

Experience flexibility

The compact photometer PF-12 is the logical result of a further development based on our very successful photometer PF-11. Adapted to our customers' requirements the PF-12 impresses with modern design and precise analytics. Pre-programmed tests, automatic wavelength adjustment and the intuitive userguidance allow fast and easy operation. Results are stored according to GLP and transferred safely to your PC with the included PC-software. The PF-12 features an unrivaled flexibility with variable power supply via mains supply, standard and rechargeable batteries, PC and 12 V car adapter. Additionally, the PF-12 offers photometric analysis of *VISOCOLOR® ECO* tests and *NANOCOLOR®* tube tests including all relevant parameters for water and waste water analysis and individual combination of tests in the new reagent cases.

Save time

Backlit graphic display with self-explanatory userguidance



- All tests and menu items can be activated fast and easily
- Operation without complex and time-consuming training

Results within seconds

The progressively designed optics is insensitive to external light and makes measuring straightforward

Be prepared for the future

Fast photometer update-free of charge

- At any time, stay up-to-date with our easy online software update
- For current software updates please visit www.mn-net.com



Assure results

Documentation of results according to GLP

- Individual entries of sample number, sample location and dilution

Clear memory management

- GLP-conform storage of results with all supplementary information such as date, time, sample number, sample location and dilution
- Fast and easy access to stored results and data sets

Convenient data export



- *NANOCOLOR®* Software DVD included in delivery
- Easy transfer of results to PC
- Data export directly to MS Excel
- Recording of calibration curves to program user-defined methods

Internal quality control according to ISO 9001

- Fulfill supervisor and authority requirements
- Fast and easy self-monitoring of photometric accuracy with *NANOCONTROL NANOCHECK* (REF 925 701)

Photometer for water analysis

Compact photometer PF-12

Be mobile

Works under any condition

- Variable power supply for mobile use:



Standard and rechargeable batteries for more than 2000 off the line measurements

- Power saving auto-off function after 5, 10, 15 or 20 min
- Waterproof housing (IP 67)



USB adaptor (REF 919 220) for cigarette lighter in the car



USB adaptor for mains supply in the lab

Enjoy versatility

Complete water and waste water analysis

- Ammonium, Chlorine, total N, total P, Nitrate, Nitrite, all COD tests and many more
- More than 100 preprogrammed tests
 - VISOCOLOR® ECO tests
 - NANOCOLOR® tube tests
- Programmable for 10 user-defined methods
- Photometric basic functions: absorbance, transmission, factor and standard



Endless opportunities

Individual mini-labs

- Portable reagent cases with photometer PF-12
- Photometric analysis of VISOCOLOR® ECO tests can be combined with:



Qualitative test papers, VISOCOLOR® reagent bottles,



QUANTOFIX® test strips, pH-Fix test strips, pH indicator papers



and accessories



Please find further information and examples of reagent cases with photometer PF-12 on page 73.

Photometer for water analysis

Compact photometer PF-12

Technical data

Type:	Filter photometer with microprocessor control, self-test and auto-calibration Wavelength range 340–860 nm
Optics:	Automatic filter wheel with 7 interference filters Insensitive to external light—no light shield required
Wavelengths:	345 / 436 / 470 / 540 / 585 / 620 / 690 nm plus 1 compartment for an additional filter
Wavelength accuracy:	± 2 nm, bandwidth at half transmission 10–12 nm
Light source:	Tungsten lamp
Detector:	Silicon photodiode
Blank value:	Automatic
Measuring modes:	Over 100 preprogrammed tests (<i>NANOCOLOR</i> ® tube tests and <i>VISOCOLOR</i> ® <i>ECO</i> tests) Absorbance, transmission, factor, standard 10 freely programmable methods
Photometric range:	± 3 E
Photometric accuracy:	± 1 %
Stability:	< 0.002 E/h
Cuvette holder:	Round tubes 16 mm OD
Data memory:	200 results, GLP conform
Display:	Backlit graphic display, 64 x 128 pixels All important data at a glance: Result in respective unit, date, time, sample number, sample location and dilution
Operation:	Display userguidance, plastic foil keyboard Test selection via test number or parameter lists 12 languages (de, en, fr, es, it, nl, hu, pl, pt, cz, id, si)
Quality Control:	With <i>NANOCONTROL</i> <i>NANOCHECK</i>
Interface:	USB 2.0
Update:	Via Internet / PC, free of charge
Operating range:	0–50 °C, up to 90 % relative humidity
Power supply:	Via USB power supply, standard or rechargeable batteries
Housing:	Waterproof, IP 67
Dimensions:	215 x 100 x 65 mm
Weight:	0.7 kg
Warranty:	2 years



This device complies with the following directives:
 - 2006/95/EG - Low-Voltage Directive
 - 2004/108/EG - EMC-Directive

Ordering information

Description	REF
Compact photometer PF-12 Incl. software DVD, manual, 4 batteries, 4 empty test tubes, funnel, beaker, syringe, USB cable, calibration cuvette and certificate in rugged case	919 200
Accu pack	919 201
USB power supply	919 220
Charger	919 221



Photometer for water analysis

Compact photometer PF-12

VISOCOLOR® ECO and NANOCOLOR® tests which can be evaluated with the PF-12

Ordering information

Test	Ranges	Test no	Wavelength	No of tests	REF
VISOCOLOR® ECO					
Ammonium 3*	0.1–2.5 mg/L NH_4^+	5-08	690	50	931 208
Ammonium 15*	0.5–8.0 mg/L NH_4^+	5-10	585	50	931 210
Chloride*	1–60 mg/L Cl^-	5-18	470	90	931 218
Chlorine 2, free + total	0.10–2.00 mg/L Cl_2	5-15	540	150	931 215
free Chlorine 2	0.10–2.00 mg/L Cl_2	5-16	540	150	931 216
Chlorine 6, free + total	0.05–6.00 mg/L Cl_2	5-17	540	200	931 217
free Chlorine 6	0.05–6.00 mg/L Cl_2	5-19	540	400	931 219
Chlorine dioxide*	0.20–3.80 mg/L ClO_2	5-21	540	150	931 221
Chromium(VI)*	0.04–1.00 mg/L CrO_4^{2-}	5-20	540	140	931 220
Copper	0.1–5.0 mg/L Cu^{2+}	5-37	585	100	931 237
Cyanide*	0.01–0.20 mg/L CN^-	5-22	585	100	931 222
Cyanuric acid*	10–100 mg/L Cya	5-23	540	100	931 223
Fluoride	0.1–2.0 mg/L F^-	5-27	585	150	931 227
Hydrazine*	0.05–0.40 mg/L N_2H_4	5-30	436	130	931 230
Iron 1*	0.10–2.00 mg/L Fe	5-25	540	200	931 225
Iron 2	0.04–2.00 mg/L Fe	5-26	540	100	931 226
Manganese*	0.1–5.0 mg/L Mn^{2+}	5-38	436	70	931 238
Nickel*	0.1–5.0 mg/L Ni^{2+}	5-40	470	150	931 240
Nitrate*	1–80 mg/L NO_3^-	5-41	436	110	931 241
Nitrite	0.02–0.50 mg/L NO_2^-	5-44	540	120	931 244
Oxygen*	1–8 mg/L O_2	5-88	540	50	931 288
pH 6.0–8.2	pH 6.0–8.2	5-70	436/540	150	931 270
Phosphate*	0.2–5.0 mg/L $\text{PO}_4\text{-P}$ 0.6–15 mg/L PO_4^{3-}	5-84	690	80	931 284
Potassium*	2–25 mg/L K^+	5-32	690	60	931 232
Silica	0.2–3.0 mg/L SiO_2	5-33	690	80	931 233
Sulfate*	20–200 mg/L SO_4^{2-}	5-92	436	100	931 292
Sulfide*	0.05–0.80 mg/L S_2^-	5-94	620	90	931 294
Zinc	0.1–3.0 mg/L Zn^{2+}	5-98	620	120	931 298
NANOCOLOR® tube tests					
Aluminum 07	0.02–0.70 mg/L Al^{3+}	0-98	540	19	985 098
Ammonium 3	0.04–2.30 mg/L $\text{NH}_4\text{-N}$ 0.05–3.00 mg/L NH_4^+	0-03	690	20	985 003
Ammonium 10	0.2–8.0 mg/L $\text{NH}_4\text{-N}$ 0.2–10 mg/L NH_4^+	0-04	690	20	985 004
Ammonium 50	1–40 mg/L $\text{NH}_4\text{-N}$ 1–50 mg/L NH_4^+	0-05	690	20	985 005
Ammonium 100	4–80 mg/L $\text{NH}_4\text{-N}$ 5–100 mg/L NH_4^+	0-08	585	20	985 008
Ammonium 200	30–160 mg/L $\text{NH}_4\text{-N}$ 40–200 mg/L NH_4^+	0-06	585	20	985 006
AOX 3*	0.1–3.0 mg/L AOX 0.01–0.30 mg/L AOX	0-07	470	20	985 007
BOD ₅ *	0.5–12.0 mg/L O_2	8-22	470	25-50	985 822
BOD ₅ -TT*	0.5–7.5 mg/L O_2	8-25	436	11-21	985 825
Cadmium 2	0.05–2.00 mg/L Cd^{2+}	0-14	540	10-19	985 014
Carbonate hardness 15	1.0–18.0 °e 0.4–5.4 mmol/L H^+	0-15	436/585	20	985 015
Chloride 50*	0.5–50.0 mg/L Cl^-	0-21	470	20	985 021
Chloride 200*	5–200 mg/L Cl^-	0-19	470	20	985 019

* This product contains harmful substances which must be specially labeled as hazardous. For detailed information please see MSDS.

visicolor®

Photometer for water analysis

Compact photometer PF-12

Test	Ranges	Test no	Wavelength	No of tests	REF
Chlorine/Ozone 2*	0.05–2.50 mg/L Cl ₂ 0.05–2.00 mg/L O ₃	0-17	540	20	985 017
Chlorine dioxide 5	0.15–5.00 mg/L ClO ₂	0-18	540	20	985 018
Chromate 5	0.05–2.00 mg/L Cr(VI) 0.1–4.0 mg/L CrO ₄ ²⁻	0-24	540	20	985 024
COD 40*	2–40 mg/L O ₂	0-27	345	20	985 027
COD 60*	5–60 mg/L O ₂	0-22	345	20	985 022
COD 160*	15–160 mg/L O ₂	0-26	436	20	985 026
COD 160 Hg-free*	15–160 mg/L O ₂	0-26	436	20	963 026
COD 300*	50–300 mg/L O ₂	0-33	436	20	985 033
COD 600*	50–600 mg/L O ₂	0-30	620	20	985 030
COD 1500*	100–1500 mg/L O ₂	0-29	620	20	985 029
COD 4000*	400–4000 mg/L O ₂	0-11	620	20	985 011
COD 10000*	1.00–10.00 g/L O ₂	0-23	620	20	985 023
COD 15000*	1.0–15.0 g/L O ₂	0-28	620	20	985 028
COD 60000*	5.0–60.0 g/L O ₂	0-12	620	20	985 012
org. Complexing agents 10 (Screeningtest)	0.5–10.0 mg/L IBiC	0-52	540	10–19	985 052
Copper 7	0.10–7.00 mg/L Cu ²⁺	0-54	585	20	985 054
Cyanide 08*	0.02–0.80 mg/L CN ⁻	0-31	585	20	985 031
DEHA 1 (Diethylhydroxylamine)	0.05–1.00 mg/L DEHA	0-35	540	20	985 035
Ethanol 1000	0.10–1.00 g/L EtOH	8-38	620	23	985 838
Fluoride 2	0.1–2.0 mg/L F ⁻	0-40	620	20	985 040
Formaldehyde 8*	0.1–8.0 mg/L HCHO	0-41	585	20	985 041
Formaldehyde 10	0.20–10.00 mg/L HCHO	0-46	436	10–19	985 046
Hardness 20	1.25–25.0 °e 5–50 mg/L Mg ²⁺ 0.2–3.6 mmol/L 10–100 mg/L Ca ²⁺	0-43	540	20	985 043
HC 300* (hydrocarbons)	0.5–5.6 mg/L HC 30–300 mg/kg HC	0-57	436	20	985 057
Iron 3*	0.10–3.00 mg/L Fe	0-37	540	20	985 037
Lead 5*	0.10–5.00 mg/L Pb ²⁺	0-09	540	20	985 009
Manganese 10*	0.1–10.0 mg/L Mn ²⁺	0-58	470	20	985 058
Methanol 15	0.2–15.0 mg/L MeOH	8-59	620	23	985 859
Molybdenum 40*	1.0–30.0 mg/L Mo (VI) 1.6–50.0 mg/L MoO ₄	0-56	345	20	985 056
Nickel 7*	0.10–7.00 mg/L Ni ²⁺	0-61	470	20	985 061
Nitrate 8*	0.30–8.00 mg/L NO ₃ -N 1.3–35.0 mg/L NO ₃ ⁻	0-65	345	20	985 065
Nitrate 50*	0.3–22.0 mg/L NO ₃ -N 2–100 mg/L NO ₃ ⁻	0-64	345	20	985 064
Nitrate 250*	4–60 mg/L NO ₃ -N 20–250 mg/L NO ₃ ⁻	0-66	345	20	985 066
Nitrite 2	0.003–0.460 mg/L NO ₂ -N 0.02–1.50 mg/L NO ₂ ⁻	0-68	540	20	985 068
Nitrite 4	0.1–4.0 mg/L NO ₂ -N 0.3–13.0 mg/L NO ₂ ⁻	0-69	540	20	985 069
total Nitrogen TN _b 22*	0.5–22.0 mg/L N	0-83	345	20	985 083
total Nitrogen TN _b 60*	3–60 mg/L N	0-92	345	20	985 092
total Nitrogen TN _b 220*	5–220 mg/L N	0-88	345	20	985 088
Organic acids 3000*	30–3000 mg/L CH ₃ COOH 0.5–50.0 mmol/L CH ₃ COOH	0-50	470	20	985 050
Oxygen 12*	0.5–12.0 mg/L O ₂	0-82	436	22	985 082
Peroxide 2	0.03–2.00 mg/L H ₂ O ₂	8-71	620	10–19	985 871

* This product contains harmful substances which must be specially labeled as hazardous. For detailed information please see MSDS.

Photometer for water analysis

Compact photometer PF-12

Test	Ranges	Test no	Wavelength	No of tests	REF
pH 6.5-8.2	pH 6.5-8.2	0-72	436/540	100	918 72
Phenolic index 5	0.2-5.0 mg/L Phenol	0-74	470	10-19	985 074
Potassium 50*	2-50 mg/L K ⁺	0-45	690	20	985 045
ortho and total Phosphate 1	0.05-1.50 mg/L P 0.2-5.0 mg/L PO ₄ ³⁻	0-76	690	19	985 076
ortho and total Phosphate 5*	0.20-5.00 mg/L P 0.5-15.0 mg/L PO ₄ ³⁻	0-81	690	19	985 081
ortho and total Phosphate 15*	0.30-15.00 mg/L P 1.0-45.0 mg/L PO ₄ ³⁻	0-80	690	19	985 080
ortho and total Phosphate 45*	5.0-50.0 mg/L P 15-150 mg/L PO ₄ ³⁻	0-55	690	19	985 055
ortho and total Phosphate 50*	10.0-50.0 mg/L P 30-150 mg/L PO ₄ ³⁻	0-79	436	19	985 079
POC 200 (polyoxycarboxylic acids)	20-200 mg/L	0-70	436	20	985 070
residual Hardness 1	0.03-1.25 °e 0.004-0.180 mmol/L	0-84	540	20	985 084
Silver 3	0.20-3.00 mg/L Ag ⁺	0-49	620	20	985 049
Starch 100	5-100 mg/L starch	0-85	540	19	985 085
Sulfate 200*	10-200 mg/L SO ₄ ²⁻	0-86	436	20	985 086
Sulfate 1000*	200-1000 mg/L SO ₄ ²⁻	0-87	436	20	985 087
Sulfide 3*	0.05-3.00 mg/L S ²⁻	0-73	620	20	985 073
Sulfite 10*	0.2-10.0 mg/L SO ₃ ²⁻	0-89	436	20	985 089
Sulfite 100*	5-100 mg/L SO ₃ ²⁻	0-90	470	19	985 090
Surfactants: Anionic surfactants 4*	0.20-4.00 mg/L MBAS	0-32	620	20	985 032
Surfactants: Cationic surfactants 4*	0.20-4.00 mg/L CTAB	0-34	620	20	985 034
Surfactants: Nonionic surfactants 15*	0.3-15.0 mg/L Triton® X-100	0-47	620	20	985 047
Thiocyanate 50*	0.5-50.0 mg/L SCN ⁻	0-91	470	20	985 091
TOC 25*	2.0-25.0 mg/L TOC	0-93	585	10	985 093
TOC 60	10-60 mg/L TOC	0-94	585	10	985 094
TOC 600*	40-600 mg/L TOC	0-99	585	10	985 099
TTC/Sludge activity 150*	5-150 µg TPF 0.050-2.300 E	8-90	470	20	985 890
Zinc 4*	0.10-4.00 mg/L Zn ²⁺	0-96	620	20	985 096

* This product contains harmful substances which must be specially labeled as hazardous. For detailed information please see MSDS.