Manufactured by MACHEREY-NAGEL
OPTIMA®
LIPODEX®
HYDRODEX
Made in Germany

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Fused silica capillary column for GC · Installation and operating instructions

Column installation	8
Start-up and testing column and chromatograph	10
General information	13
Technical support and customer services	14

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Visit our chromatography pages

www.mn-net.com

- Free application database with more than 3000 chromatography applications
- Derivatization protocols
- Integrated internet shop
- Product information
- Material safety data sheets
- Job opportunities

The screenshot displays the MACHERY-NAGEL website's chromatography section. At the top, there is a navigation bar with links for Home, SPC, HPLC, TLC, GC, Flash/Absorbents, Vials/Closures, Filters, Service, Applications, and WebShop. Below this, a 'Products' section features icons and descriptions for various categories: SPE (CHROMABOND® HR-XenT concept, SPE accessories), HPLC (NUCLEODUR®, NUCLEOHELL®, NUCLEODUR™ Special phases), TLC (Laves' Introductory kits, Absorbents), GC (OPTIMA® PERMABOND®, Dimensional kits, Special phases), Flash LC absorbents (Flash chromatography, Irregular silica, LC absorbents), Vials and caps, CHROMAFIL® filters, and Applications. A 'Service' section offers Support, Technical customer support, and a link to the website. The main content area is titled 'MN Application database HPLC - GC - TLC - SPE' and states 'Now more than 3000 applications available!'. It includes a search interface with fields for 'Search criteria: Substance / Keyword', 'Search word', and 'Application type' (with checkboxes for SPE, TLC, HPLC, GC). A 'search' button is located below these fields. On the right side of the search area, there is a small chromatogram showing several peaks. The footer of the website includes the copyright notice 'Copyright 2013 by MACHERY-NAGEL GmbH & Co. KG' and links for 'Home', 'Impressum', 'Anwendungsdienste', and 'Kontaktierungsmöglichkeiten'.

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Please read these instructions BEFORE installing – and certainly before heating – your column!

The present capillary column from MACHEREY-NAGEL has been manufactured in our own facilities and will be the heart of your gas chromatographic analytical system. It has been produced following state-of-the-art technologies and has been individually tested; the enclosed certificate documents its quality standard.

In order to obtain good results from your column as long as possible, you should carefully read and follow these instructions. We have grouped them into 4 chapters:

- I. Column installation
- II. Start-up and testing column and chromatograph
- III. General information for handling
- IV. Technical support and customer services

I. Column installation

The following items are required for proper column installation:

- ◆ Ferrules (with hole appropriate for the respective column outer diameter)

Column ID [mm]	Column OD = bore of ferrule [mm]	REF (graphite ferrule)	Notation of other manufacturers
0.05–0.20	0.40	708309	minibore
0.25	0.40	708309	narrowbore
0.32	0.50	708308	widebore
0.53	0.80	708301	megabore

For GC instruments of Carlo Erba / Fisons and Hewlett-Packard / Agilent special ferrules are available – please see our website www.mn-net.com or our catalogue Chromatography.

- ◆ Diamond file (e.g., REF 708300)
- ◆ Magnifying lens (e.g., REF 706296)
- ◆ Vial with volatile solvent (e.g., pentane, cyclohexane, acetone)
- ◆ Leak detector or leak check liquid
- ◆ Flow meter

The great number of our customers and their very different state of chromatographic knowledge make it impossible to write instructions, which cover all demands. Thus we have tried to strike a balance between an extensive textbook and succinct keywords and hope that this will help the beginner as well as the experienced user.

Important note: Do not remove the cords from the capillary column! The cords attach the column to the cage and prevent unnecessary contact between the column and metal parts of the cage. The cords are temperature stable up to 400 °C.

Now you are ready for installation:

Caution: Avoid sharp bends. The column should be installed in an unstrained configuration. A column hanger should only contact the cage. For coupling of columns use only zero dead volume fittings or devices to avoid peak broadening, tailing or even memory effects.

- a) Open the plastic packing (cut – not tear) and remove the column.
- b) Check the column for mechanical transport damage; there must be no free ends visible (ends not melted)! If this should be the case, return the column to your distributor for replacement.
- c) The melted column ends have to be removed. For this purpose, score the polyimide layer about 5 cm from the melted end using a diamond file and break the end off straight. Make sure that no column fragments get into the column. For this purpose it is advisable to take off the removed end in a downward movement.
- d) Unwind both column ends from the cage for about one turn.
- e) Suspend the column in the middle of the oven.

- f) Place a column nut and a ferrule over each column end, then cut about 5 cm from each column end using a diamond file (in order to remove any ferrule particles from the column).
- g) Examine the column ends with a magnifying lens. If the column ends are not smooth or even show strong burrs, cut off another piece from the column.
- h) Consult the manual of your GC instrumentation because it is very important to install the column entrance at the correct insertion distance in the injector. Insert one column end into the injector and attach it according to the instructions of the instrument manufacturer, i.e., fix it only slightly more than finger-tight.
- i) Adjust the flow of the carrier gas through the injector according to the instrument manual, the analytical instructions or the average values of the following table:

ID	Length	Pressure	Approximate dead time
[mm]	[m]	[bar]	[min]
0.05	10	2.5 (H ₂)	1.5
	20	4.5 (H ₂)	3
0.20	12	0.8 (He)	1.0
	25	1.5 (He)	1.8
	50	3.0 (H ₂)	3.0
0.25	10	0.3 (He)	1.0
	15	0.4 (He)	1.2
	25	0.7 (He)	2.5
	30	0.8 (He)	2.8
	50	1.5 (He)	5
	60	1.8 (He)	5.5
0.32	10	0.3 (He)	0.8
	15	0.4 (He)	0.9
	25	0.5 (He)	1.4
	30	0.6 (He)	1.6
	50	1.0 (He)	3.5
	60	1.2 (He)	4.0
0.53	10	0.2 (He)	0.6
	15	0.2 (He)	0.7
	25	0.3 (He)	1.0
	30	0.4 (He)	1.1
	50	1.0 (H ₂)	1.5

- j) Place the free end of the capillary column into a vial with a volatile solvent. After not more than 30 s gas bubbles should indicate a steady flow through the column. If no bubbles appear, repeat steps f) – i). If you still don't see a gas flow check the following:
- the supply of carrier gas
 - the pressure of the carrier gas in the supply pipes (pressure regulator, valves)
 - if the capillary column is broken
- k) Insert the second column end into the detector and attach it according to the instructions of the instrument manufacturer, i.e., fix the nut slightly. Exception: If using highly sensitive detectors (MS, ECD) you should first condition the column (see below) and then connect it to the detector.

Hint: For subsequent column installations you may mark the installation depth of the column in the injector and detector in the gas chromatograph using a moisture-proof felt tip pen (e.g., on the door or the edge of the oven).

- l) Check the column connections in the oven for leaks (e.g., with leak check liquid; this is especially important when using hydrogen as carrier gas)!
- m) Close the oven door, start the gas chromatograph:
- Heat the injector and the detector (maximum to the long-term temperature limit of the capillary column) and ignite, if required (e.g., FID).
 - Heat the oven (choose a temperature between 60 and 100 °C).
 - Start the integrator or the data processing system.
 - Adjust the split ratio (typical: 1:30).
 - Inject methane or a volatile solvent to check the overall performance of the gas-chromatographic system.
- n) According to requirements, condition the column (over night), (see chap. II.1).
- o) Possibly reproduce the column test from certificate (see chap. II.2) – now the column is ready for use.

II. Start-up and testing column and chromatograph

1. Conditioning

Each column has been conditioned prior to its quality control. Conditioning is used to elute volatile components of sample molecules and stationary phase and to remove active sites of the inner column surface.

Conditioning should always be performed after a longer storage of the column.

Conditioning can be performed by repeated running of a temperature program (3–5 cycles) or a longer exposure (12–24 h) to the maximum isothermal temperature.

The two numerical temperature values mentioned in the certificate indicate the (lower) isothermal long-term temperature (= conditioning temperature) and the (higher) short-term temperature (= maximum in temperature gradients during 5–15 min).

During conditioning never exceed the long-term temperature, and make sure that the gas supply is never interrupted.

Caution: *Without carrier gas flow at high temperatures any stationary phase in a capillary column will be destroyed in no time!*

2. Column test

Each column from MACHERY-NAGEL is tested in a gas chromatograph, and the individual chromatogram with conditions (the test certificate) is supplied with the column. It should be retained for future reference.

If possible you should reproduce the included test chromatogram before using the column for your analytical tasks. The respective test mixture can be separately purchased (see www.mn-net.com).

This will indicate whether

- ◆ column installation is correct
- ◆ samples are properly injected
- ◆ the overall chromatographic system including evaluation works accurately.

By reproducing the column test you document the performance of the column, and thus you can at any time compare the state of the column with the initial quality.

Using the column test, you can identify a number of typical failures; this is the reason why reproduction of this test should always be the first step in troubleshooting.

If the test result is satisfactory, the planned application can be operated.

3. Checking column installation

What should you do when your own test does not show the same result as the test certificate?

- a) Check the correct installation of the column and test again.
- b) Compare the faulty chromatogram with the examples in our troubleshooting and try to correct the problem.
- c) Contact your nearest MACHERY-NAGEL office or your local distributor and ask for help.

In the following chapter we have listed some typical shortcomings which can be seen in the chromatogram.




4. Troubleshooting




If you cannot reproduce the test chromatogram of your capillary column, this may be due to one of the following reasons.

For example, a dead volume in the injector produces tailing. This is especially strong for early eluting components in a temperature program. Dead volume in the detector also results in tailing, however, in this case all components are affected.

The following "faulty" chromatograms are by-products of years of laboratory work; suggested remedies are based on experience with a large number of problems from daily practice. We have selected failures which from our point of view are most frequently encountered, but of course our choice can by no means be comprehensive.

Fused silica capillary column for GC

Observation	Possible cause	Suggested remedy
Shortcoming A 	Sample too diluted Column connection leaks, column not properly installed Injector temperature too low Sample decomposes in the injector Column adsorbs substances or decomposes them	Increase injection volume, concentrate sample Check for leaks, change ferrules, reinstall column Check temperature program, increase injector temperature Check temperature program, reduce injector temperature, change liner, check capillary ends Check capillary ends, check intact deactivation using a suitable test mixture; for poor results shorten both column ends by about 10 cm; or replace column If the column test does not show any defects: a) use a column with thicker film; b) use phase with better deactivation; c) use special column
Missing or overlapping peaks		
Shortcoming B 	Septum particles in the column entrance Heating program too fast, final temperature too high Column bleeding due to incomplete column conditioning Soiled detector	Cut off one turn from the column entrance, change injection septa Decrease temperature gradient and final temperature Disconnect detector and condition column (also see instructions for your GC instrumentation and detector) Clean detector according to manufacturers instructions
Increasing baseline, at high temperatures bleeding, strong noise		
Shortcoming C 	Bleeding of silicone septa	Use new septum or septum grade with lower bleeding (high temperature septum)
Regular interfering peaks		
Shortcoming D 	Impurities from sample vials and / or sample preparation Dirty syringe Sample is decomposed Column adsorbs substances or decomposes them	Check SPE and autosampler vials, use septa with lower bleeding Use clean syringe Check temperature program, oven temperature (thermometer); when chromatographing thermolabile compounds, reduce injector temperature, change liner Check capillary ends, check intact deactivation using a suitable test mixture; for poor results shorten both column ends by about 10 cm; or replace column If the column test does not show any defects: a) use a column with thicker film; b) use phase with better deactivation; c) use special column
Irregular interfering peaks / spikes / ghost peaks	Bleeding of butyl septa	Use silicone septa for injector and sample vials

<i>Observation</i>	<i>Possible cause</i>	<i>Suggested remedy</i>
<p>Shortcoming E</p> 	<p>Column overload Sample is decomposed</p> <p>Column adsorbs substances or decomposes them</p>	<p>Reduce injection volume</p> <p>Check temperature program, oven temperature (thermometer); when chromatographing thermolabile compounds, reduce injector temperature, change liner</p> <p>Check capillary ends, check intact deactivation using a suitable test mixture; for poor results shorten both column ends by about 10 cm; or replace column</p> <p>If the column test does not show any defects: a) use a column with thicker film; b) use phase with better deactivation; c) use special column</p>
<p>Fronting, heavy peak broadening of the peak front</p>	<p>Sample with high boiling point System leaks Sample is decomposed</p> <p>Column adsorbs substances or decomposes them</p>	<p>Derivatize polar, basic or high-boiling samples</p> <p>Check for proper column installation and leaks, change ferrules</p> <p>Check temperature program, oven temperature (thermometer); when chromatographing thermolabile compounds, reduce injector temperature, change liner</p> <p>Check capillary ends, check intact deactivation using a suitable test mixture; for poor results shorten both column ends by about 10 cm; or replace column</p> <p>If the column test does not show any defects: a) use a column with thicker film; b) use phase with better deactivation; c) use special column</p>
<p>Shortcoming F</p> 	<p>Compounds which always show strong tailing</p>	<p>No remedy possible</p>
<p>Tailing, heavy peak broadening of the peak back</p>	<p>Compounds which always show strong tailing</p>	<p>No remedy possible</p>
<p>Shortcoming G</p> 	<p>Split flow too low Column overload</p>	<p>Increase split flow</p> <p>Inject less or increase split flow</p>
<p>Broad peaks</p>		

III. General information for handling

1. Column regeneration

If contamination is the cause of a loss in column performance (compare shortcoming D) you may try one of the following remedies:

- a) Remove the contaminated end from the column, i.e. in general shorten the injector column end by about 50 cm (at least one column turn).
- b) Remove column end in the transfer line of the detector, because here phase deterioration may be faster than in the rest of the column due to the permanent high temperatures.
- c) You may also try to "flood" the column in the GC using 20 μ L acetone followed by a normal heating program.
- d) As last resource the column can be rinsed (*only chemically bonded phases*). For this purpose disconnect the capillary column from the gas chromatograph and connect the detector end (*mark it!*) to the rinsing device. For rinsing use 10 mL pentane, acetone or toluene for nonpolar – medium polar phases or dichloromethane for polar phases. Never use acids or bases (*even diluted*) or derivatization reagents! After rinsing dry the stationary phase for at least 2 h in a stream of carrier gas at ambient temperature before you install the column in the GC and heat it in a temperature program (1 $^{\circ}$ C/min) up to the maximum isothermal temperature. Hold this temperature for about 1 h.

2. Column storage

If you want the capillary column to stay in the GC over night we recommend to reduce the split flow and adjust the oven temperature to 80–100 $^{\circ}$ C.

If you have to turn off the carrier gas or if it is necessary to change the column, the oven should be cooled to ambient temperature before disconnecting the column. In order to prevent oxygen damaging of the stationary phase the column should be stored with sealed column ends (*especially important for Wax and FFAP columns*). Fused silica columns can be easily sealed by sticking the ends (sideways) into a septum. For long-term storage we recommend flame sealing and storage in the original shipping container.

3. Maintenance of value

The column at hand consists of a tube of fused silica manufactured with high precision and coated with a protective polyimide layer which results in a high column flexibility.

Caution: *It is very important to avoid damage to the polyimide layer, since unprotected fused silica becomes brittle and thus easily breaks.*

The column is coiled onto a cage which largely protects it from mechanical damage. Avoiding mechanical or extreme thermal strain (e.g., opening the oven door at high temperatures!).

In addition, you can protect your column from chemical stress by performing a suitable sample preparation, by controlling the pH value of your sample (neither acid nor basic) and by using a precolumn. If matrix problems occur, a change in sample preparation can considerably extend the lifetime of the column.

We recommend that you use CHROMAFIL[®] syringe filters for sample clarification, CHROMABOND[®] SPE columns for selective sample preparation and/or deactivated capillary columns (precolumns) to protect the separation column from contamination. If necessary 0.5–1 m (always at least one column turn) of the guard column can be removed without need to shorten the analytical column.

IV. Technical support and customer services

If you have any questions concerning the use of capillary columns please don't hesitate to contact us:

MACHEREY-NAGEL GmbH & Co. KG

Germany and international	Phone	+49 24 21 969-0
	Fax	+49 24 21 969-199
	E-mail	info@mn-net.com
Switzerland	Phone	+41 62 388 55 00
	Fax	+41 62 388 55 05
	E-mail	sales-ch@mn-net.com
France	Phone	+33 388 68 22 68
	Fax	+33 388 51 76 88
	E-mail	sales-fr@mn-net.com
USA	Phone	+1 484 821 0984
	Fax	+1 484 821 1272
	E-mail	sales-us@mn-net.com

We will be glad to help you with:

- Technical support by phone
- Technical documentation as comprehensive collections of applications for our products (e.g., GC Application Guide)
- Local technical support (e.g., by your local distributor)
- Our extensive internet collection of applications for GC, HPLC, TLC, SPE
- Visit our chromatography pages: www.mn-net.com

Liability

MACHEREY-NAGEL products are intended for general laboratory use only. They are not suited for clinical use (diagnostic, therapeutic or prognostic). MACHEREY-NAGEL does not assume any responsibility for damages due to improper application of our products or due to application of our products in other fields of application. The user has to ensure that the products used are suitable for the planned application. Application on the human body is strictly forbidden! The user is liable of all damages resulting from such application.

MACHEREY-NAGEL products are suited for qualified personnel only.

All goods have to be checked immediately upon receipt. Damage claims are only acceptable in writing within 8 working days of receipt of all the goods. In case of legitimate claims the buyer can only require replacement of the goods. If replacement is not possible, the buyer has the right to choose between alternative products with same value or refund. The buyer cannot claim further compensation. All returns must first be authorized by us in writing. We offer a one-year warranty from date of delivery that our products will conform to applicable specifications set forth in the product specifications if not sold to persons set forth in § 13 of German Civil Code. In such a case, the provisions of the German Civil Code shall be valid.

This warranty does not cover defects in or damage to the products which are due to improper installation or maintenance, misuse, neglect or any use other than ordinary commercial application.

Any discharge from liability will be void if a defect results from a heavily negligent or intentional breach of contract or if the buyer will be bodily injured as a result from a negligent or intentional breach.

The seller shall not be liable for the goods being fit for a particular purpose unless otherwise agreed upon, to which the buyer intends to put them. This warranty is strictly exclusive. Any further damage compensation is impossible.

MACHEREY-NAGEL makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER; DIRECTLY OR INDIRECTLY; EXPRESS OR IMPLIED; INCLUDING; WITHOUT LIMITATION; AS TO THE SUITABILITY; REPRODUCIBILITY; DURABILITY; FITNESS FOR A PARTICULAR PURPOSE OR USE; MERCHANTABILITY; CONDITION; OR ANY OTHER MATTER WITH RESPECT TO MACHEREY-NAGEL PRODUCTS OR THE SALE OF MACHEREY-NAGEL PRODUCTS.

If not superseded by these terms the terms and conditions for the sale of goods (export version) are additionally valid.

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MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany
MACHEREY-NAGEL AG
MACHEREY-NAGEL Inc.