**Elektrische Spindelpressen – Motor-driven pumps
Pompes électriques**

**300–7000 bar
4300–100 000 PSI**

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
**Model A
Motor-driven pumps with adjustable flow**

- For the generation of pulsation-free and constant flow.
- Controlled constant volumetric flows also at minimum flow (control range 1–100 %).
- Prolonged constant flow due to large stroke capacity.
- Hand operation (pulse key forward and reverse) or automatic, continuous discharge flow.
- Quick return to starting position and automatic reverse-to-forward motion.
- Pressure generator with asynchronous positioning motor with integrated frequency converter and encoder (2048 pulses/revolution).
- Combination of inductive and mechanical final position switches for reversing of movement and stroke limitation.
- Control box 500×400×210mm (h×w×d) with 3.5" TFT touch panel (indication of actual volumetric flow and position). Supply 230 VAC/50 Hz.
- Digital interface for start/stop and direction of movement and analogue interface for set point (volumetric flow: standard configuration 0–10V).
- Analogue output for actual volumetric flow and piston position (standard configuration +/-10V).
- All wetted parts in stainless steel and PTFE.

Options

- Different analogue signal for external set point or actual value (e.g. 0/4–20mA).
- Model B: With pressure control.
- Model C: Without controls.
- Opposed-piston pump (boxer).

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
- Equipment like model A, but with pressure control function.
- To maintain a constant pressure in a high-pressure system.
- Included equipment: pressure controller, pressure transducer, hand valves and check valve (according to flow diagram).
- Digital indication of piston position.
- Digital interface for start/stop and direction of movement and analogue interface for set point (pressure: standard configuration 0–10 V).
- Analogue output for actual pressure (standard configuration 0–10 V).

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Model B Motor-driven pumps with pressure control

- Design like model A but with pressure control function.
- To maintain a constant pressure in a high-pressure system.
- Included equipment: pressure controller, pressure transducer, hand valves and check valve (according to flow diagram).
- Digital indication of piston position.
- Digital interface for start/stop and direction of movement and analogue interface for set point (pressure: standard configuration 0–10 V).
- Analogue output for actual pressure (standard configuration 0–10 V).

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

Alternating

Continuous discharging of liquids. While the first piston is discharging, the cylinder volume of the directly coupled second piston is filled. It will thus be ready for its following discharge stroke.

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**Model C
Motor-driven pumps without controls**

- Assembled on a base plate with reduction gear box, asynchronous positioning motor and limit switches (an inductive and a mechanical limit switch on each end position).
- Asynchronous positioning motor with integrated frequency converter and encoder (2048 pulses/revolution).
- Serial programming interface and corresponding software is in the scope of delivery.
- Programmable digital and analogue inputs.

Option

- Motor-mounting adapter according to the requirement of the customer.

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Compensating

For measurements of solubilities gas/liquid and liquid/liquid with sampling. In a test cylinder the discharge volume of the first piston compensates at constant pressure the suction volume of the directly coupled second piston.

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