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Instructions for use

For correct test results, carefully follow the instructions for use given below. Do not perform the test at temperatures below +15° and above 35°C.

Test procedure

Follow this sequence exactly:

1. Keep watch with seconds hand ready. Note down date and time.
2. Open container. Remove one test strip. Do not touch the reaction zone!
3. Reseal container immediately after use.

Test strips for rapid determination of glucose in blood

Application

Self-monitoring of glucose concentrations in capillary blood for diabetics resp. suspected diabetics. For monitoring of diabetics evaluation systems with higher accuracy should be applied. **As far as diabetes mellitus is concerned, this test in no way replaces regular professional medical controls.**

Drawing blood sample

The blood sample is obtained by pricking the fingertip. First wash hands with warm water to improve blood circulation and facilitate sampling. **Do not use the first drop of blood.** Please read the instructions for use of your lancet system.

Measuring range

20–800 mg/dL (1.1–44.4 mmol/L)

Results

If evaluation of the test field 60 seconds after wiping the blood sample (120 seconds after application) yields a value greater than 120 mg/dL, the glucose concentration in the blood is too high. Consult your doctor.

Never change the intake of prescribed medication on the basis of the results obtained with these test strips.

Directions

Discuss any **positive** (glucose concentration >120 mg/dL) or any **doubtful** test results with your doctor. This, obviously, is also required if test results are negative, but pain persists and/or you do not feel well.

If a falsely negative test result is suspected, repeat the test using a test strip from an as yet unopened container or consult your doctor. Your doctor will then conduct further examinations using the methods available in his/her laboratory. In any case, in order to establish a final diagnosis and prescribe the appropriate therapy, your doctor should verify results obtained with test strips, with diagnostic methods of his own.

Remove only as many test strips as required. Do not touch the test fields! Reseal the vial immediately after use.

Do not remove the dessicant inside the cap, do not damage the cap.

Protect the test strips from sunlight and moisture. Store the vial in a cool, dry place (not in the refrigerator) between 4 ... and 30°C. When stored correctly the test strips can be used up until the expiry date printed on the vial, also after initial opening.

Throw away damaged test strips or vials.

The caps contain a non-poisonous and harmless dessicant. Should this dessicant be swallowed, then drink plenty of water (consult a doctor if necessary).

Dispose of used test stripes (test stripes can only be used once) along with household waste.

Explanation of symbols on folding box and container label see below. Keep test strips out of reach of children!

Explanation of symbols

 Statement of Conformity (Product corresponds to the In-Vitro Diagnostic Medical Devices Directive 98/79/EC of the European Union)

 Please read instructions for use!

 Permitted storage temperature range

Use by

 Batch identification

 Item number

 Manufacturer

 Do not reuse

Information and additional notes**Principle**

Detection is based on the glucose oxidase-peroxidase reaction and reacts specifically to glucose.

Evaluation – Sources of Error

Visual evaluation by comparison with the color scale. Normal concentrations (fasting blood sugar) are between 70 and 120 mg/dL (3.9–6.7 mmol/L). (Comparison values measured on Ebio plus (Eppendorf))

Self-test studies have shown that 95 % of all readings have been assigned to the correct color field. 3 % of the readings showed a falsely positive result, only 2 % of all readings were falsely negative. Both hyperglycemia and hypoglycemia can be determined with this test. If the test results vary slightly from normal values, a quantitative (e.g. photometric) test should be carried out. Glycaemie C is **unsuitable** for monitoring neonatal blood specimen.

A weaker color change is to be expected after intravenous injection of metformin and ascorbic acid and with administration of dopamine. The effect of medicines or their metabolic products on the test is not known in all cases. In case of doubt, consult your doctor.

Extreme hematocrit values (< 30 and > 55 %) influence the result. Severely lipemic (fatty) blood samples interfere with the detection reaction: in case of fat metabolism disorders, please consult your doctor. When venous blood samples are drawn, only EDTA shall be used as anticoagulants. Sodium fluoride (NaF) causes falsely negative results.

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